

Product datasheet for **MC212448**

Lrrc3c (NM_001195544) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lrrc3c (NM_001195544) Mouse Untagged Clone
Tag: Tag Free
Symbol: Lrrc3c
Synonyms: Gm12
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212448 representing NM_001195544
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCCTACTAGCTCCTCGCCTCCTGCCCTCCTGCTGGTGATTGGCACCGGGGCTCTGTGCCCCGCC
 CCCTGGCACTGCCCCAGGGCTGCTACGTAGCAGAAGAAGCCGGGAACAGACATTCGGCTGCAGCCGGGC
 TGGTCTGAGTGCCGTGCCAACGGTATCCCAACGACACCCGAAAGCTTTACCTGGACGCCAACAGTTG
 GCATCAGTGCCAGCTGGTGCCTTCCAGCACCTGCCTGCCCTGGAAGAAGTACCTGTCTCATAACGCTC
 TTGTTCTCTCTCAGGGGCTGCTTCCAGGGCCTGGAGGGCACTCTGCGCCACCTCGACCTCTTGCCAA
 CCAGCTAGCATCTGTGCCTGTGGCAGCCTTCGTGGGGTGCAGATCCAAGTGAACCTGTCTGCCAACCCA
 TGGCGCTGCGACTGTGCCTACAGGAAGTGCTCAGGCACGTGAGTTGGCTCCAGGCTCCGGGACAGGCA
 TTGTGTGGTCCAGAAGCCCGACCAGATCTCGTGGGACACGAGTTCTTCTGCTCACCAGGGAGGAAGA
 GCTGTGTGGGACAGGACGAGGTGGGACCCGAAGGAGCACGGATGTGGCCCTGTTAGTCACCATGGGGGC
 TGGCTGACACTGGTAGTGGCTTATCTGATCCGCTATGTGTGGCAGAACCAGATGAGACGCGCGCCCGG
 TCAAGCGGGCTCCTCCGCGCAGCCTGTGCGCTCAGAGGACTCCTTACACTCAGCACAGTGCTC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001195544
Insert Size: 768 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001195544.1</u> , <u>NP_001182473.1</u>
RefSeq Size:	902 bp
RefSeq ORF:	768 bp
Locus ID:	192852
Cytogenetics:	11 D