

Product datasheet for MC212450

Cyb5d2 (NM_001024926) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cyb5d2 (NM_001024926) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cyb5d2
Synonyms:	9330151E16Rik; Gm2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC212450 representing NM_001024926 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCTGAGGATATGCGGGCTTGGTGTGGTGTGAGCCTGGCCGTGGCGCGGTGGCTGTCATGGCCGTGT
GGTTGATGGACTGGTGGGGTCCCCGCCCTGGTATTCGCCTCTTCCTTCCCGAGGAGCTAGCCGCTATCG
CGGCGGCCAGGGGACCGGGCTGTACTTGGCACTGCTCGGCCGCGTCTACGACGTGCTCTTGGCAGG
AGGCACTACGAGCCTGGGGCGCACTATAGCGGCTTCGAGGCCGGGATGCATCCAGAGCATTGTGACAG
GAGACTATTCTGAAGCCGGTCTTGTGGATGATATAAATGGCTTGCCTCCTCTGAGATACTGACGCTGCA
CAACTGGCTTTCACTTCTACGAGAAAAATTATGATTTGTCGGAAGACTGGTTGGAAGTTTACAGAAAG
GATGGGTTACCCACTTCAGAACTAACCCAGGTAGAAGCCATGGTGAATAAGGCATGGAGGCAAAATGAAC
AGGAACAAAGAGAAAAGCAGAAGTTCCACCATGCAATTCGAGTGGAGCTCTGCTAAGGGCAGCAGGCT
CTGGTGTCTCCAAAAGAGTGGAGGTGTGCACAGAGACTGGATTGGTGTCCCCAGGAAGTTGTATAAGCCA
GGGGCCAAGGAGCCCCATTGTGTGTGTGTGAGAACAACCTGGCCCTCTAGTGACCAGCAAGACAACCCTA
GACACTCAAATCATGGGGACTTGGACAACCCCACTTGAAGAATACACAGGCTGCCCACTGGCTAC
CACATGTTCTTCCCACTCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001024926
Insert Size:	792 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001024926.3, NP_001020097.2</u>
RefSeq Size:	2393 bp
RefSeq ORF:	792 bp
Locus ID:	192986
UniProt ID:	<u>Q5SSH8</u>
Cytogenetics:	11 B4
Gene Summary:	Heme-binding protein which promotes neuronal but not astrocyte differentiation. [UniProtKB/Swiss-Prot Function]