

Product datasheet for MC212473

Zfp691 (NM_183140) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp691 (NM_183140) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Zfp691
Synonyms:	6430544H17Rik; Znf691
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC212473 representing NM_183140 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGCAGCGAGAAGGAGCAGCGTCCCGAAGCTCACCTGCCCGAGGAAGGAGAAGGGGCTAAACCTTGGAGAGTGGATGGCTCAAAGGATTCTCAGATCACACCTCGGGAGGATCATGGGCAGGAGAGCCTGTTGGCAGGCTCCACGGAACGCATCCACCAAAGACAAGGCAGAAAGTCACTGCCCAAGCCGGAGGCCCCAGGGATCCCATGCTTTTTTCAAGCCCAGAGACAGATGAGAAGCTTTTTATATGTGCGCAGTGTGGCAAAACCTTCAACAATACCTCAACCTGAGAACGCACCAGCGGATCCACACTGGCGAGAAGCCCTACAAGTGTTCGAGTGTGGCAAGAGTTTCTCCGGAGCTCCAACCGCATCCGGCAGCAGCGCATCCACCTGGAAGAGAAGCACTACCAAAGCAGTGTGTCAGGAGAGCTTCCGGCGCGCTCGGACCTCACCACGCACCAGCAAGATCACCTGGGCCAGCGGCCCTACCGCTGTGACATTTGTGGCAAGAGCTTCACTCAGAGCTCCACGCTGGCTGTACATCACCGGACCCACCTGGAGCCAGCCCCTTACATCTGCTGTGAGTGTGGGAAGAGCTTTAGCAACAGCTCCAGCTTTGGGGTGCACCACCGTACCCACACGGGTGAGAGGCCTTACGAGTGCAGTGAATGCGGGCGGACTTTCAGTGATATCTCCAACCTTGGTGCACACCAGAGAACGCACAGAGGGGAGAAGCCGTACAGGTGCACCCTGTGTGGAAACTTCTCCGGAGCTCCAACCTCATCCGACACCAGAAAACACACTTGGCGAGCAAGATGAGAAGGATTTTAGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_183140
Insert Size:	852 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_183140.3 , NP_898963.2
RefSeq Size:	2338 bp
RefSeq ORF:	852 bp
Locus ID:	195522
UniProt ID:	Q3TDE8
Cytogenetics:	4 D2.1
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript. All three variants encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.