

Product datasheet for **MC212478**

Cetn4 (NM_145825) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cetn4 (NM_145825) Mouse Untagged Clone
Tag: Tag Free
Symbol: Cetn4
Synonyms: MGC107235
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212478 representing NM_145825
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCATCCAGCCAGCGCATAACTTTAGACCAATGGAAGAAAAAGCAGCAAAAGTTGAACTGAATGACA
CCCAGAAGCAAGAAATTAAGAGGCCTTTGACTTATTCGATATTGATGGATCTGGAACCATAGATCTGAA
AGAACTGAAGATTGCAATGAGGGCCTTAGGCTTTGAGCCAAAGAAAGAAGAAGTAAACAACCTGATCGCT
GAAATCGACAAGGAAGGGACTGGACCATTGTTTTGAAGATTTTTTGCATCATGAGCGTAAAAATGA
GTGAAAAAGATGAAAAAGAAATACTGAAGGCATTTAAATTATTTGATGACGATGCTACGGGAAGCAT
ATCACTAAATAACATCAAGAGGGTTGCCAAGGAGCTAGGGGAAAACCTAACCGAAGATGAATTACAGGAA
ATGCTTGATGAAGCTGATCGTGACGGGGATGGAGAGATTAATGAGGAAGAATTCTTGAAGATGATGAAGA
AGACCAGCCTTTAT**AG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_145825
Insert Size: 507 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145825.2 , NP_665824.1
RefSeq Size:	1059 bp
RefSeq ORF:	507 bp
Locus ID:	207175
UniProt ID:	Q8K4K1
Cytogenetics:	3 B
Gene Summary:	Ca(2+)-binding protein that may be involved in basal body assembly or in a subsequent step of ciliogenesis.[UniProtKB/Swiss-Prot Function]