

Product datasheet for MC212517

Zc3hav1l (NM_172467) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zc3hav1l (NM_172467) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Zc3hav1l
Synonyms:	B130055L09Rik; E430016P22Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC212517 representing NM_172467 Red=Cloning site Blue=ORF

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGACCTCACTCTGTGCGGTTCTCACCAAGGTGCTGTGTGCTCACGGCGCCGCATGTTCTGCG
AGGACCTGCGCGCCACGTGGAGCTGTCGGAGGCGAAGCTGCGGGCCGTGCTGCGCGCGCGGGTCCAGA
GCGCTTCTGCTGCAGGAGGTGGAGCTGCGCGACGGCCCGTGGGACGAGAGGCCGAGGTGGCGCGGGC
GAGGGCGCGGGCGGGCTCAGGCGCGGCCACGGCTGCAGGGTGATGGCCGTCTCTCCGCGCGCC
TCTGCGCCGCTACCAGCGCGGGAGTGTGAGGGCTGCGACAGCTGCACTTGTGTGCGGCCACATGCT
GGCAAGTGCCCCACCGGACTGCTGGTCTACCTGTTCTGTGCCATGATATCCACACACCTATCAAC
GTCCAAGTCTGAAAAGCCGGGGCTTTTTGGCCTCAATGAAGGGCAGCTGCGAATCCTGCTTTTGAGA
ACGACCCCTGCCTTTCCAGAGGTCTGTGAGATGTACAACAAAGGTGTCGATGTCCTTTATGGCTACTG
CAGCCTCAAGGATAGATGCAACAAGTTTCACGTGTGCAATCCTTTGTGCGAGGTGAATGCCCTTTTCAG
CCCTGCAACGCTCCACAGCTCATCCATGCTGCCACCTCAAGTTGCTAGAGGACCAAGAACTTAGTG
TTTCAAGCGTCGTGAATTTCCAGATAATCTGTCTACAGACATAAGAAGTGCACAAGATGCTTGAAGA
GAAAGACCACTCGGCTTCTACTGAGCAACCACAAGGCCTTGGGAAACAAGGAGCCCTCGGGGCTGTGAA
GCCAGGCTTTCTCCCTGCCCGTGTCTAGTCGCCAGGAAGCCCAATAA

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms:	https://cdn.origene.com/chromatograms/ja1719_g09.zip
Restriction Sites:	SgfI-MluI
ACCN:	NM_172467



[View online »](#)

Insert Size:	891 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	BC138180 , AAI38181
RefSeq Size:	1981 bp
RefSeq ORF:	891 bp
Locus ID:	209032
UniProt ID:	Q8BFR1
Cytogenetics:	6 B1