

Product datasheet for MC212591

P2ry10b (NM_175442) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	P2ry10b (NM_175442) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	P2ry10b
Synonyms:	A630033H20Rik; AI662791; P2ry10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC212591 representing NM_175442 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGTAGCAACAGTACCAGCAGTGCTGAGAGTAATTGCAATGGCACTTATTTGACATTTCAATACTCTC
TGTATGCAACCACCTATATCTTCATATTCATCCCGGTCTCTGGCTAACAGTGTAGCCCTGTGGTTCT
GTGCCGTTTCATCAGCAAGAAGAACAAGGCCATCATCTTCATGATCAACCTCTCAGTGGCAGACCTTGCT
CATGTCCTGTCTTACCTCTTCGGATTTACTATTATATCAATCGTCACTGGCCATTCCAGAGGGCCCTTT
GCCTGCTGTGTTCTATCTGAAATATCTCAACATGTATGCCAGCATTTTTTCTTGACATGCATTAGCCT
TCAGAGGTGCCTTTTCTCTCAAGCCATTTCAGAGCCAGAACTGGAAGCGTAGGTATGATGTGGCTATC
AGTGCTGTCTGTGGATCATTGTGGGACTGCCTGTTTGCCACTCCCATCCTGAGAAGTCTGGCTTAG
CCAATAACAGTGAATCGTGCTTTGCTGATTTAGGACTTCAGGACATTACTATGGCTTCATCCATTGGCAT
GGTAACTGCAGCTGAACCTGGAGGGTTGTATTGCCTGTTGTAATTATTACATATTGCACATGAAAAACA
AGAAAATCTTTACAGGAATTCAGATCCACCTCAAAATATTAAGAGAGGAAAAAGGCTTTGAGAATGG
TTCTAATGTGTGCAGTGGTATTCATAGTCTGTTTCACTCCTTACCATCTCAACTCCCATTCTTTATGAT
GGTAAAGCAACGTGTCTTTCAAACGTTCCTTTATTAGGAACACTCTATGTTCCACATCATTTTCCTTG
TGCTTTGCAAATCTGAATTGTTGCTTGATCCAGTTGTGTTATTTATGACTTCAGAAATTCGTGATC
AATTTTCAGAACATGGGAATTTGGTCAATCCAGTCATGTTTGCATTGTAAGGACAGTACTTTGGAATTCG
TCAGAGGAAGGAGAATCTTCAGACTATCTCTGGAATGTTGGATGTTCCAACAATGTGATGAAATG
GTTATCTAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_175442
Insert Size:	1059 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_175442.4, NP_780651.1</u>
RefSeq Size:	3008 bp
RefSeq ORF:	1059 bp
Locus ID:	213438
Cytogenetics:	X D