

## Product datasheet for **MC212661**

### Gm11545 (NM\_001105561) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Gm11545 (NM_001105561) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gm11545
Synonyms:	A430060F13Rik; C86086; OTTMUSG00000002043
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC212661 representing NM_001105561 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCTGGACACCTGGGTCCCGGGCAGGCTCACCTTGATTTTTGTCCTGATTTCCAGTCTCCGACGAGTGA  
GCACCAACGAGACTGTAACACATGTGACATCTTAACTGCCCAAAGGATTCACATGTTGTGACAAGGA  
CTGCTGCCGAAAGAAAGGTCTGGGGTCTGAAAGCTGGAATCCTAAAAATGAACCTTCAAGATCGTG  
GCATCATCCTATTGGTCATGCTACCCCTTTATGCATCTGCGGTCTAGTGAGGCCTTCTGTCCCAATT  
GTAGAGAACTACAGCATAACTTCAGAACAGCAGATCATCAGACATCTCCATCGTTGCCCTCCGTTGCTCC  
CCTGGAGAACATATGGGTACCTCCTTGGATCCCCACCAGCCTACAGTCAGGTTGTTCTGAAGCCAAC  
CCCACAGAACCACCTCCTCCCTACAGCCTCAGGCCTGAGGGCACTTCTGGACAGAGGAGAGGGCATGCCT  
ACGCAACCCTT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001105561
Insert Size:	504 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001105561.1, NP_001099031.1</u>
<b>RefSeq Size:</b>	2707 bp
<b>RefSeq ORF:</b>	504 bp
<b>Locus ID:</b>	217122
<b>Cytogenetics:</b>	11 D