

Product datasheet for MC212711

Apol9a (NM_173786) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Apol9a (NM_173786) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Apol9a
Synonyms:	MGC68363
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC212711 representing NM_173786 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAAGCAGTGGATACCTGTATAGCCTGGTGGCTGAAGATGGAGTCTGGCAGGCGTTTGTATCCAAGG
CCAAGTTGTCCAGGGAGAGAGCAGTGGCCCTCCACAAAGCTCTGAGGGAAGTACAGCGCTCTTGGCCAT
AGCAGACAGAGGCAGACTTCAAAAAGGTCTGAAAGGCAGGGAGAAGTTTTGAAAGTCTTTCCTCGCTTG
AAAGCAGAGCTGGTGGAGCATATCAGTCAGCTCCATGCCCTGGCTGACCACGCTGAGAACTACACAGG
GCTGCACCATCTCAACGTGGTGGCTGACTCCTCAGTGCTGCCTCTGACATCCTGAGCCTCCTTGGTCT
CTTCTGGCACCTGTGACAGCAGAGGGAAGTCTGGTGCTCTCAGCAACTGGCTTGGGGCTGGGGGTAGCA
GCTACTGTGACTGATGTTGCTACCTCAATCGTGGAGGAAACAAGCAGGTTTTGGATGAAGTTGAAGCTG
CTCCGACTGGCACTCATGTGCTGGGGGAGGCTGGCACAAGTGTGGATAGGATTGCCAGCAAGGTCCTCA
GGCTACCAGAGACATCACTAGAGACCTGGAAGCCCTTGAGCAGCACATGAATGCCCTCAGGCTGGTTAGA
GCCAACCTCGTTTAGAAGAAGATGCCAGGATCCTGGCCACCACAGGAAGCATCTCTGCCAACGGGCTA
TGCAGGTGCGGGCCAGTCTGAAAGGAACCCCTCTGGCAATGAGCAATGAAGCCGGATCCAACGTGCCAC
CGCCGCAGGTGCCGCCCTCTGGAGTGATGTGGACAGCCTTGTCAAAGAGTCAAAGCATCTGTACGAAGGG
TCTGCGTCAGAGTCGGCTGAAGCACTGAGGAAGCTGGCTCGAGAGCTGGAGGAGAAGCTAGAGGGGCTCA
TGAATTCTACAAGACGATC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_173786



[View online »](#)

Insert Size:	933 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_173786.1</u> , <u>NP_776147.1</u>
RefSeq Size:	1527 bp
RefSeq ORF:	933 bp
Locus ID:	223672
Cytogenetics:	15 D3