

Product datasheet for MC212717

Gtsf2 (NM_177631) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gtsf2 (NM_177631) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gtsf2
Synonyms: MGC58362
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212717 representing NM_177631
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGTTGGAATCCATAGAGCTATGTCCTTATGATCCAAACCACAGGATACCTGCCAGCAGGTTACAGT
ACCATCTGGCATCATGCAAAAAGAAAAATCAAAGATAGCTAAAAAGATGGCTAACTGCAATATAATGC
TTGTCATGTGGTCCCAATCAAAGGCTCAAGGAACATGAGGCCAACTGTATCAACAGGACTGCAGTGGAT
GATGAGCCACTTAACCTTCAAAGATTACTCACCCAGTTTTTGAAGAAAATGAAAACCTTCCAGTGGCT
GGAGTCAATTTCTGACCTGATGTCTGGAACGTGGATCATGCGCATCATTTTCCTTCATTTGTTCTTGA
GACTTTTGCTCCAAAGAAGCTAGTTTGTGAGAGTGACTCGAGAGACCTACAAGAGGCTATGGCTGATAAG
CATCCTAACAGCTTCAAGTCTTGGGGAAGAGGGCAGAAAAC TGA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_177631

Insert Size: 465 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_177631.3, NP_808299.1</u>
RefSeq Size:	753 bp
RefSeq ORF:	465 bp
Locus ID:	223927
Cytogenetics:	15 F3