

Product datasheet for **MC212724**

E330017A01Rik (NM_175011) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: E330017A01Rik (NM_175011) Mouse Untagged Clone
Tag: Tag Free
Symbol: E330017A01Rik
Synonyms: MGC117520
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212724 representing NM_175011
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTTTTCTGCCATCCCTTCTCTGAAGAATGCACAAAGGCTTTGAACCAAGTTGTGGCTTACCACCTGC
ACACCTCCCATGTGTATTTTGAATGGCCACAGTTTCATGACTGGCAAGAAGGAGACTTTTCCTTTTGT
CCAGCACTTTGAAACCCTGGCTGATTCAAGGAGGAATATGCAAATAGATTCTAAACCATCTATGGACC
CGTAACAAAAAGATCTGCCCTCCTCAACGTATGAGAAGTTGATATGAAGGAAATAACTACACCTAAAA
GTGCCCTACTAATGGCCAGAATATGGAAGAGACATTAACCAAAATCCTGCTAGCCCTAAAGGCTGCTGC
ACGTCGTGAATCAGATCTCTAAAATTTTTAACATCGGTCCTCTGTAACAAAGAGAGTTTAAAGAGTAT
CTGGAAAGTCACTTTCTCCTCCGCAAAAAGGTAAAAGAAAGAATGAGAAGGAAGCCAGACTGAACAAC
CTTCTACTTCATCGGGTGTAAAAGGTACCTCAGTAGCAAAAAGAAGCAAGCCTGGAGATCATAAAGCCAC
CACCACCTCCAATTG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_175011
Insert Size: 579 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_175011.2, NP_778176.1</u>
RefSeq Size:	677 bp
RefSeq ORF:	579 bp
Locus ID:	224247
Cytogenetics:	16 C1.2