

## Product datasheet for **MC212730**

### **D17Wsu92e (NM\_001033279) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** D17Wsu92e (NM\_001033279) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** D17Wsu92e  
**Synonyms:** AU020189; C030010A15; C80239  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >MC212730 representing NM\_001033279  
**Red=Cloning site Blue=ORF**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGAGGGCATGGACGTGGACTTGGATCCGGAGCTCATGCAGAAGTTCAGCTGCCTGGGCACCACAGACA  
AGGACGTGCTCATCTCCGAGTTCACGGCTACTCGGCTTCCAGCTCAACCCGGCCGGCTGCGCCTTCTT  
CCTGGACATGACCACTGGAATTACAAGCAGCGATTGGCGCTATTATGACTTCGAGAGCCAAATATC  
AGTGTGCCCTCCATGTCTTTGTTGAAGATGTACCATAGGAGAAGGGGAGTCAATACCCCTGATACCC  
AGTTCATAAAAACTTGGCGGATACAGAATTCTGGAGCAGAGGCCCTGGCCTCCAGGGTTTGTCTTAATA  
TGTCCGGGGAGACCAGTTTGGACATGTGAACATGGTGATGGTGAGATCACTAGAGCCCAAGAGATCGCA  
GACGTCAGCGTCCAGATGTGCAGTCCCAGCAGAGCAGGCATGTACCAGGGACAGTGGCGGATGTGCACTG  
CTACTGGACTCTACTATGGAGATGTCTCTGGGTGATTCTCAGTGTGGAGGTGGGTGGACTTTTAGGAGT  
AACGCAGCAGCTGTATCTTTTGAACGGAATTAACACACAGCCACATCGCAAGGTAGAAGGAACTTC  
AACCCGTTTGCCTCTCCCAAAAGAACCACAATCAGATGAAAACAACCTAACAGACCCTGGGGTTCCG  
AGTTCGACTCGATCAGCAAAAATACATGGGCTCCTGTTCTGAGCAAAGCGAGCAAGATCAGGATAGACT  
GTCACAGAGCTGTGAAATCTGTCCCCAGCAGTCCCGCAAACAATTATCAGTAGTGACTTACAGTAAG  
GGGCTCCACGGGCTTACCCCTTCGGCCAGTCTTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001033279  
**Insert Size:** 876 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">BC145403</a> , <a href="#">AAI45404</a>
<b>RefSeq Size:</b>	1169 bp
<b>RefSeq ORF:</b>	876 bp
<b>Locus ID:</b>	224647
<b>UniProt ID:</b>	<a href="#">Q3TT38</a>
<b>Cytogenetics:</b>	17 14.59 cM
<b>Gene Summary:</b>	<p>Negative regulator of innate antiviral response. Blocks IRF3-dependent cytokine production such as IFNA, IFNB and TNF. Interacts with IRF3 and inhibits IRF3 recruitment to type I IFN promoter sequences while also reducing nuclear levels of the coactivators EP300 and CREBBP.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>