

Product datasheet for **MC212804**

Lcn11 (NM_001100455) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lcn11 (NM_001100455) Mouse Untagged Clone
Tag: Tag Free
Symbol: Lcn11
Synonyms: Gm109
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC212804 representing NM_001100455
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCAAGAAGACACAGCAGAGGAATGAAGCTTCTACTGCTGCTGAGTGTGGGGCTGGCCTGACCTGGG
 CCCTGCAGGACTCCACCCTGAGCAGGTAACAGGACCTTGGCACACCCTCAAGCTGGCTTCCACTGACCG
 GTCTCTCATTGAGGAAGGAGGTGCCTATCGGTGCTTCATGACTGACATTGTCTCCTGGACAATGGCAAC
 CTGAATGTACCTACTTCCACAGGAAAGATGAAAAGTGTGTGAAGGAATTCTACATCGCAGAGAAGACTG
 ACACACCTGGACAGTACACCTTTGAGTACCAAGGCAGAACTCTCTGACTTTTGTGCATGTTACTGAGGA
 CTTGCGCCATCATGGACTTGGAAAACAGAGTGAAGGGGGCACCACCATTGTAATTGAGTTCCACGGCAGG
 AGTCTGTCCACAGATGAGCTCGGATGGGAACGCTATCTCGTTCACACACGCAGAAGAGGCATCGCACCA
 AGAACATTGTGGACCTTTCCTTATCTCGTCGATGTGACACACACTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001100455
Insert Size: 537 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100455.2, NP_001093925.1</u>
RefSeq Size:	652 bp
RefSeq ORF:	537 bp
Locus ID:	227630
Cytogenetics:	2 A3