

## Product datasheet for MC212813

### Olfr358 (NM\_207235) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Olfr358 (NM\_207235) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Olfr358  
**Synonyms:** GA\_x6K02T2NLDC-33697309-33696317; MOR159-4  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC212813 representing NM\_207235  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGCCACTAAAAACAAGACAGAAGTGACTGAGTTTGTGTTGTTGGGATTGTCCAGTAGGCCAGAGATGC  
 AGCCAGTTATTTTTGGAGTTGTCCTTATCATGTACCTGGCAGCAGTTTTGGGCAATACCTTACTAGTGCT  
 TGTTGCTTGCTCAGACCCAGGCTTCAGACCCCATGTACTTCTCTCAGCCAGCTTTCCTTGATTGAC  
 ATATGTCTAACAACCATCACTGTTCCCTCAGATGCTGGTACACACACTCTCTGTGGTTCTGAAGCATTCTT  
 ATAAGTGTTCATGACACAATTATTTTTCTTTATGGCAGTGGGCAGCATGGAAGGCCACTTACTGGCTGC  
 CATGGCATATGACCGCTATGTTGCCATCTGTGACCCTTAAAGATACTCTGCTATTGTGAGCCATAGCCTC  
 TGTCTGCGAATAACATTGACCTCATGGTGGTGGTCAGCCTAACAGTCTCCTGTACAGTGTGTTGGTCA  
 CTCGTTAACTTCTGTGGCAACCAGGTCACCCACTTCTCTGTGACATCACTCCTCTGTGAAGCTCTC  
 TTGACTCGACCAAGTGGTCAATGAAATGTTAATTTTCACTGAGGGTGTAGCTGTGGTTGTGAGCCCTTT  
 TTCTTCATTTGGGCTCCTACATTGCGATAGGTTTTGTTATTGCCACATGCATTCAACTGCTGCCCTGA  
 GCAAAGCTGTCTACCTGTGGCTCCACATTATGGTTGTGTTGCTCCTGTATGGTTGATTGCGAT  
 GTACCTCAAGCCATCGTCTACCTATGACTTAAACCAGGATCGACAGATTGCCATCTCTACACATTGGTC  
 ACCCCTATGCTAAATCCATTAATCTATAGCCTCAGGAACCAGGAAGTTAAAGGAGCTCTCACACGGCTTC  
 TTAGGAAACTCTGCATCTCAGGAAGCTTCCAGGTTGGCTCAGATAAAGAAAGAATGACAAATATTATGAG  
 ACCTAATGTT**TAA**

**ACGGT**ACGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_207235



[View online »](#)

<b>Insert Size:</b>	993 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_207235.1</a></u> , <u><a href="#">NP_997118.1</a></u>
<b>RefSeq Size:</b>	993 bp
<b>RefSeq ORF:</b>	993 bp
<b>Locus ID:</b>	227789
<b>Cytogenetics:</b>	2 B
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]