

## Product datasheet for **MC212835**

### Tgif2 (NM\_173396) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tgif2 (NM\_173396) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tgif2  
**Synonyms:** 4921501K24; 5730599O09Rik; C80753; C81206  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC212835 representing NM\_173396  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGTCGGACAGCGATCTAGGCGAGGATGAGGGGCTCCTGTCCTGACTGACCGCAAGAGGAAGCGCAGAGGGA  
ACCTGCCCAAGGAGTCAGTAAAGATCCTCCGGGACTGGCTCTATCTGCACCGCTACAACGCCTACCCCTC  
AGAGCAGGAGAAGCTAAGTCTCTCTGGACAGACCAACCTCTCGGTGCTGCAGATATGTAAGTGGTTCATC  
AATGCCCGTCGTCGCCTTCTCCCTGACATGCTTCGGAAGGATGGCAAAGACCCTAATCAGTTCACCATCT  
CCCGTCGTGGGGTAAGGCCTCTGATGTGGCTCTCCCCAGGGGCAGCAGCCCTCTCTGCTGGCTGTGTC  
TGTTCCAGCCCCACTAACATGCTCTCCTGTCTGTGTGCTCCATGCCACTCCACTCAGGCCAAGGTGAA  
AAGCCAGCAGCACCCCTTCCCACAAGTGGAGCTGGAATCCCCAAGGCCCTGGTGACCCCTGCTAGCACAC  
TCACCCTGCTGACTAGGGCAGAGGCTGGAAGCCCCACGGGTGGACTCTTCAATACGCCACCACCCACACC  
CCCAGAGCAGGACAAGGATGACTTCAGCAGCTTCCAGTTGCTGGTGGAGGTGGCGCTGCAGAGGGCTGCC  
GAGATGGAGCTTCAGAAGCAGCAAGAGCCAGCGCCACCTTGTACACACTCCGCTGCCTTTCGTCTCAG  
AAAACGCCAAG**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_173396  
**Insert Size:** 714 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_173396.3</a> , <a href="#">NP_775572.1</a>
<b>RefSeq Size:</b>	2900 bp
<b>RefSeq ORF:</b>	714 bp
<b>Locus ID:</b>	228839
<b>UniProt ID:</b>	<a href="#">Q8C0Y1</a>
<b>Cytogenetics:</b>	2 H1
<b>Gene Summary:</b>	<p>Transcriptional repressor, which probably repress transcription by binding directly the 5'-CTGTCAA-3' DNA sequence or by interacting with TGF-beta activated SMAD proteins. Probably represses transcription via the recruitment of histone deacetylase proteins (By similarity). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (a).</p>