

## Product datasheet for **MC212840**

### Stx16 (NM\_001102423) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Stx16 (NM\_001102423) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Stx16  
**Synonyms:** 4930401D03; 5430410K23Rik; 6330500A18Rik; AI648908; AW553605; Syn16  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC212840 representing NM\_001102423  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCCACCAGGCGTTTAAACCGACGCTTCTTGTGTGGCGGAATAATCCATCCAACCCGGCAGCTGT  
TAGCCGAGCAAGTAGTAGTACACACCCTCCAGCCCTCTGCATTACGTAGCATTGCTGCGCTTGTCTGA  
TGACCGCATGGCGCTGGTATCAGGCATCAGCCTAGATCCAGAAGCAGCAATTGGCGTGACAAAGCGGTCCG  
CCTCCAAAGTGGGTGGATGGCGTGGACGAAATACAGTATGATGTTGGCCGATTAAACAGAAGATGAAGG  
AGTTAGCCAGCCTTCATGATAAGCATTTGAACAGACCCACCCCTGGATGACAGCAGTGAGGAGGAACATGC  
CATCGAGATCACCAACAGGAAGTACACAGCTCTTCCATAGGTGCCAGCGTGTGTGCAGGCCCTGCC  
AGCCGGGCACGAAGGGCCTGCTCGGAGCAGGAAGAGCGGCTACTGCGGAACGTGGTGGCCTCCTTGGCAC  
AGGCCCTGCAGGAGCTGTCCACCAGCTTCCGGCATGCACAGTCTGACTACCTGAAACGCATGAAAAACCG  
AGAGGAGAGATCGCAGCATTTCTCGACACCCAGTACCACTAATGGATGATGGAGATGATGCTACTCTG  
TATGGTCAGGGCTTACGGATGACCAGCTGGTGTGGTAGAGCAAAACACTGATGGTGGAGGAGAGGG  
AGCGAGAGATCCGTGAGATTGTGCAGTCCATTTAGACCTCAATGAAATCTTCAGGGACTTGGGAGCCAT  
GATTGTGGAGCAGGGAACAGTCTTGTAGAAATGACTATAATGTTGAACAGTCCCTGTGTCAAGACCGAA  
GATGGCTTGAACAGCTTACAAGGCTGAGCAGTACCAGAAGAAGAACCAGGAGATGCTTGTGATCGTGA  
TCCTGGTGGCCGTATCATAGTCTCTCGTGGCCCTCGTCGCCGTAAGTCTCGCTAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001102423  
**Insert Size:** 969 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001102423.1</a></u> , <u><a href="#">NP_001095893.1</a></u>
<b>RefSeq Size:</b>	3756 bp
<b>RefSeq ORF:</b>	969 bp
<b>Locus ID:</b>	228960
<b>Cytogenetics:</b>	2 H4
<b>Gene Summary:</b>	SNARE involved in vesicular transport from the late endosomes to the trans-Golgi network. [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (b) has the same N- and C-termini but is shorter compared to isoform a.