

Product datasheet for **MC212880**

Ifna13 (NM_177347) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ifna13 (NM_177347) Mouse Untagged Clone
Tag: Tag Free
Symbol: Ifna13
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC212880 representing NM_177347
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTAGGCCTTGTGCTTTCCTGATGGTTTTGGTGGTGTGAGCTACTGGTCAGCCTGCTCTCTAGGAT
GTGACCTGCCTCAGACTCATAACCTCAGGAACAAGAGAGCCTTGACTCTCCTGGAACAAATGAGGAGACT
CTCCCTCTTTCCTGTCTGAAGGACAGAAAGGACTTTGGATTCCACAGGAGAAGGTGGATGCCCAGCAG
ATCAAGAAGGCTCAAGCCATCCCTTTTGTTCATGAGCTGACCCAGCAGATCCTGACCCTTTCACATCAA
ATGACTCATCTGCTGCATGGAATGCAACCCTCCTAGACTCATTCTGCAATGACCTCCATCAGCAGCTCAA
TGACCTCAAAGCTTGTCTGATGCAGCAGGTTGGGGTGCAGGAATTTCCCTGACCCAGGAAGACTCCCTG
CTGGCTGTGAGGAAATACTTCCACAGCATCACTGTGTACCTGAGAGAGAAGAAACACAGCCCTGTGCC
GGGAGGTGGTCAGAGCAGAAGTCCAGAGAACCTGTCTTCTCAGCCAATTGCTGGCAAGACTGAGCAA
GGAGGAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_177347
Insert Size: 570 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_177347.2</u> , <u>NP_796321.1</u>
RefSeq Size:	819 bp
RefSeq ORF:	570 bp
Locus ID:	230396
UniProt ID:	<u>Q80SU4</u>
Cytogenetics:	4 42.02 cM
Gene Summary:	Exhibits antiviral activity against Theiler's virus, Mengo virus and vesicular stomatitis virus. Interferons alpha stimulate the production of two enzymes: a protein kinase and an oligoadenylate synthetase.[UniProtKB/Swiss-Prot Function]