

## Product datasheet for **MC212907**

### **Fbxo44 (NM\_173401) Mouse Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fbxo44 (NM\_173401) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Fbxo44  
**Synonyms:** 5730411K09; AV001623; FBG3; Fbx6a; FBX30; Fbxo6a  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC212907 representing NM\_173401  
**Red**=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGCAGTAGGCAACATCAACGAGCTGCCAGAGAACATCCTGCTGGAAGTGTTCATCCACATCCCAGCTC  
GCCAAGTCTGCTGCGCTGCCGACCGGTCTGCAGCCTCTGGCGAGACCTCATTGACCTGGTTACACTCTG  
GAAGCGCAAGTGCCTTCAAGAGGGCTTCATCACCGAGGACTGGGACCAGCCTGTGGCTGACTGGAAGATC  
TTCTACTTCTGCGCAGCCTCCAGAGGAACCTCCTTACAAACCTTGTGCTGAAGAGGGATTTGAGTTCT  
GGAGCCTGGACGTGAACGGAGGGATGAATGGAAGGTAGAGGATCTCTCTAAAGACCAGCGAAAGGAATT  
CCCCAATGACCAGGTCAAGAAATACTTCGTGACTTCTATTACACCTGCCTCAAGTCCCAGGTGGTGGAC  
CTCAAGGCTGAAGGGTATTGGGAGGAGCTGATGGACACCACCGGCCGGACATCGAGGTCAAGGACTGGT  
TTGCAGCCAGGCCAGACTGCGGGTCCAAGTACCAGCTGTGCGTCCAACCTCTGTCGTGAGCCACGCACC  
ACTGGGAACCTTCCAGCCGACCGGTGATGATTACAGCAGAAAAGCGATGCCAAGTGGAGGGAGGTTTCA  
CACACATTTCCAACATCCTCCCGCGTCCGCTACATCTGGTTTCAGCACGGAGGCGTGGACACTCACT  
ACTGGGCCGGCTGGTACGGCCCGAGAGTCACAAACAGCAGCATCACCATCGGGCCCCCGCTGCC**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_173401  
**Insert Size:** 768 bp



<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_173401.5</a> , <a href="#">NP_775577.2</a>
<b>RefSeq Size:</b>	1863 bp
<b>RefSeq ORF:</b>	768 bp
<b>Locus ID:</b>	230903
<b>UniProt ID:</b>	<a href="#">Q8BK26</a>
<b>Cytogenetics:</b>	4 E2
<b>Gene Summary:</b>	Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (1). Both variants 1 and 3 encode the same isoform.