

## Product datasheet for MC212956

### Vmn1r178 (NM\_206868) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids

**Tag:** Tag Free

**Symbol:** Vmn1r178

**Synonyms:** Gm473; V1rd13

**Mammalian Cell Selection:** Neomycin

**Vector:** pCMV6-Entry (PS100001)

**E. coli Selection:** Kanamycin (25 ug/mL)

**Fully Sequenced ORF:** >MC212956 representing NM\_206868  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGCTCTGCTCATGATAAATCTCTGAAAACCACTGATGACATGACTCTTCAAATCCTCCTGCTTTGCCAGT  
TTGTGGTTGGGACTGTGGCCAATGTCTTTCTGTTCTCCATAATATCTCTCCAGTCTTGACTGGTTCTAA  
ACAGAGGCCAGACAGTTGATTTTAAGCCACATGTCTGTGGCCAATGCCTTGATTCTCTTCTCACTTCA  
TTTCCAGACATGATGGCTTTTGCTCCAAGAAATCTCCAATCACCTCAAATGTAAATTGGAATACTTTA  
GTCGCCTGGTGGCAAGAAGCACAAACATGTGTTCCACCTGTGCTCCTGAGTATCCATCAGTTTGCACTCT  
TGTTCTGTGAATAGAGGTAAAGCTATACTCAGAGCAAGTGTCCAACTTGACGAGTTATTCTTGTAT  
AGTTGTTGGTTCTTCAGTGTCTTAAGTAATATCTACATTGCAATTAAGGTCACAGGTCCACAGATTACAG  
ACAATAACACTGACTCTAAAAGCAAGTTGTTCTGTTTCACTTCTGGATTCACTGTCGGCATTCTTCTCT  
GCGGATTTTCCATGATGGCACATTCATGAGCATCATGGTCTGGACCAAGTGTCTCCATGGTACTTCTCCTC  
CATAGACATCGCCAGCGAATGCAGCACATCCTCACTCCCAATCAGGATGCCAGAGGCCAAGCTGAGACCA  
GAGCAACCCATACTATCCTGATGCTGGTGGTCACATTTGTTAGCTTTTATCTTCTAAATTTTATCTGTAT  
CATCTTTTACACCTTTTCTATAGATTCTCGTCTCTTCATGAGTCATTTTCAGTAAAGTTTTTGCTGCAGGA  
TTCCCCACTATCTCCCTTACTGTTGATCTTCAGAGATCCTAAGTATCCTTGTCTGTGCTCTTCAACT  
GTTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_206868



|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Insert Size:</b>           | 915 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                          |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                    |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol> |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>RefSeq:</b>                | <u><a href="#">NM_206868.1</a></u> , <u><a href="#">NP_996751.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>RefSeq Size:</b>           | 915 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>RefSeq ORF:</b>            | 915 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Locus ID:</b>              | 232959                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Cytogenetics:</b>          | 7 A3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |