

Product datasheet for **MC213031**

Klhl15 (NM_001039060) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Klhl15 (NM_001039060) Mouse Untagged Clone
Tag: Tag Free
Symbol: Klhl15
Synonyms: 6330500C13Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC213031 representing NM_001039060
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCAGGGGACGTGGAAGGATTCTGTTCTCCATCCATGACACCAGTGTCTCTGCTGGGTTCCAGAGCAC
TGTATGAGGAGGGATTGCTTCTCGATGTCACCTGTTTATTGAAGACCATCAGTTCCAGGCCATAAAGC
TCTCTTGGCCACCCAGAGTGATTACTTCAAGATTATGTTACAGCAGACATGAGGGAACGAGATCAGGAC
AAAATTCATTTGAAAGGTCTAACAGCTACTGGTTTCAGCCACGTTCTTCAGTTTATGTATTATGGAAC
TAGAACTAAGTATGAACACTGTACATGAAATCCTTCAGGCAGCAATGTATGTTTCAGCTTATAGAAGTGGT
GAAATCTGCTGCTCCTTCTCTTAGCCAAAATCTGCTTAGAAAATTGTGCGGAAATATGAGACTGTTA
GATGATTTCCGGTGTAAACATCGAGGGAGTCAGGGAGAAGCTGGATGCCTTCTGCTAGACAACTTCGTGC
CACTAATGTCCAGGCCTGACTTTTTGTCTTATCTGAGCTTCGAGAAGCTCATGCTTATTTGGATAATGA
CCATCTGAGCAGGTTTCCAGAAATTGAGCTGTACGAGGCTGTGCAGTCTTGGCTGCGACATGATAGACGA
CGGTGGAGACACAGATACCATCATTCAGAACATCAGGTTTTGTCTGATGACTCCATCCAGCGTCTTTG
AGAAGTTGGTACATTTGAAAGCGATAGACAGGACTCCTCAATT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001039060
Insert Size: 747 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001039060.1</u> , <u>NP_001034149.1</u>
RefSeq Size:	4001 bp
RefSeq ORF:	747 bp
Locus ID:	236904
UniProt ID:	<u>A2AAX3</u>
Cytogenetics:	X C3
Gene Summary:	Substrate-specific adapter for CUL3 E3 ubiquitin-protein ligase complex. Acts as an adapter for CUL3 to target the serine/threonine-protein phosphatase 2A (PP2A) subunit PPP2R5B for ubiquitination and subsequent proteasomal degradation, thus promoting exchange with other regulatory subunits and regulating PP2A holoenzyme composition. Acts as an adapter for CUL3 to target the DNA-end resection factor RBBP8/CtIP for ubiquitination and subsequent proteasomal degradation. Through the regulation of RBBP8/CtIP protein turnover, plays a key role in DNA damage response, favoring DNA double-strand repair through error-prone non-homologous end joining (NHEJ) over error-free, RBBP8-mediated homologous recombination (HR).[UniProtKB/Swiss-Prot Function]