

Product datasheet for MC213100

Fam69c (NM_173770) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Fam69c

Synonyms: B230399E16Rik

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC213100 representing NM_173770
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGTGGCCGGGAGGTCAAGAACACGCTAGGCTTGGAGCTGCCCAACAACAGCATAGCCCCGCTATGGC
CTGCCAGGCAGGGCCCTGGCTGGCGGCAACAGCTGGCCAGTGCCTGCTCCAAACAGGAGGAATA
CGTATACTTCAGCCTGCTTCCAGATCTGAGCCGCCACATCCTACCGGTGCTGGGATCCTGCGGTCACTTT
TATGCAGTGGAGTACTTGGCTGCTGGTAGCCCTCACCACAAGGCTCTCTCCCACTGGATGATGCTGGCC
AAGCACAGGCTATTAGCCACATTGCCCTCAGCTTCTGGACATGGTGAGCCATTCGACAGTGACTTTTC
CCACCGGCTCCACCTCTGTGATGTCAAACCTGAGAACTTCGCCATCAAGAGGGACTTTACGGTGGTAGCT
ATTGACGTAGACATGGCCTTTTTTGGAGCCTAAATGAGGGAAATCCTTGAGCAGAACTGCACGGGAGATG
AAGATTGCAATTTCTCGATTGTTTCTCAAGTGTGATCTCCGAGTGCACAAGTGTGGAGCCAGCGTGT
CAACAGCAACCTGCAAGTCATCTGTGATAAAATATTCCGACACTGGTTCTCCTCAACACACAGGAGCCCG
GCTGTCTCCCTCCAGTTCGGCTGCAGTACAGCAGGCCGTGCAGGAGTGTGCTCAGCACGGGGGCTCCT
CGGGGAATTCCTGGACAGCCTCATCCAGTGTGTCTGGAAGCTTCGATGGCTCTTGACAGGCCACTCTGAA
AGAACTGCAGGAGCGGAGAAGTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_173770

Insert Size: 795 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_173770.4, NP_776131.1</u>
RefSeq Size:	2153 bp
RefSeq ORF:	795 bp
Locus ID:	240479
UniProt ID:	<u>Q8BQT2</u>
Cytogenetics:	18 E4