

Product datasheet for **MC213164**

Oas1f (NM_145153) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Oas1f (NM_145153) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Oas1f
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC213164 representing NM_145153 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGAAGGATCTTAGCAGCACCCAGCCTGTGAGCTGGACAAGTTCATACGTGATCATCTCCTCCCCG
ATCCAGCTTCCATGCTGAGGCCAGAGCAGACGTGGACTTCATAGGTGCTTTCCTGAAGGAGAGATGCTT
CCAAGTGCCACCCACCCTGTGAGGGTCTCCAGGTTGTGATGGGCGGCTCCTACGACGAACACACTGCA
CTCAAGAGCAAGTCAGAGGCTAAAATGGTGGTGTCTTAACAATCTCACCAGCTTCGAGGAGCAGTTAA
AGCGACGGGGAGAGTTTCATTGAGGAAATTCGGAACACCTGTGTGAGCTGCAGGATGAGAAACCATTAA
AGTGAAGTTTGGGTGCAGAGCTCAGAGGAGCCCAACTCCAGGTCTCTGAGCTTCAAGCTGAGCTCCCT
GAGCTCCAGCAGGAGGTGGAATTTGATGTGCAGCCAGCCTATGATGTCCTGTATGAACTGAGAAACAACA
AGTATGCTGAACCTACTTGTACAACAAAATCTACGCCCAACTCATCCATGAGTGCACCACACTAAAGAA
GGAGGGCGAGTTCTCCATCTGCTTACCAGCCTCCATCAGAGCTTCTGGAGGATCGTGCACCCAAGCTG
AAGAACCTCATCCGTTTGGTCAAGCACTGGTATCAACTGTGAAGGAGAAGCTGGGAAGCCGCTGCCCC
CACAGTATGCCCTGGAGCTGCTCACAGTGTACGCCTGGGAAAGTGGGAGTAGAGACTGCGAATTCACAC
AGCCCAGGGCTTCGAACTGTCTTGAACTGGTACCAAGTACAAGTGGCTTCGAATCTACTGGACAGTG
TATTATGACTTTAGAAAGACGAAGTCTCTGAATACCTGCACAAAATGCTCCAAAAGTCAAGCCTGTGA
TCCTGGACCCTGCTGACCCAACAAGGAACGTGGCTGGTACCAACCTACTAGGCTGGGGGCTGTTGGCAAA
AGAAGCTGCCATCTGGCTGCAGTCTCCTGCTTTAGGAACTGTGATACGTGCCTCGTGGGCCCTGGGGT
GTGCCGATGAAGGTCGAGATCCACAGGACTGTGTCCTTCTA**TGA**

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_145153
Insert Size:	1095 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145153.3 , NP_660135.2
RefSeq Size:	1480 bp
RefSeq ORF:	1095 bp
Locus ID:	243262
Cytogenetics:	5 60.64 cM