

## Product datasheet for MC213165

### Hcar1 (NM\_175520) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Hcar1 (NM_175520) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Hcar1
Synonyms:	Gpr81
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC213165 representing NM_175520 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCCAGTCTCTCCAAGTCTATGGACAACGGGTCGTGTCTCATCGAGGGGAGCCCATCTCCC  
AAGTGATGCCTCCTACTCATCCTGGTCTTCGTGCTTGGCGCCCTGGGCAACGGCATAGCCCTGTGCGG  
CTTCTGCTTTCACATGAAGACCTGGAAGTCAAGCACTATTTACCTTTCAACTGGCTGTGGCCGATTTT  
CTCCTCATGATCTGCTTACCCCTTCGGACAGACTACTACCTCAGACGCAGACACTGGATTTTTGGAGATA  
TCGCCTGTGCGCTGGTCCTTCAAGCTGGCCATGAATAGGGCCGGGAGCATTGTCTTCTCACTGTGGT  
GGCTGTGGATAGGTATTTCAAAGTGGTCCACCCACCATATGGTGAATGCCATCTCCAACCGGACTGCC  
GCCGCCACCGCTGTGCTCTGGACTTTGGTCATCTTGGGACTGTGTATCTTCTGATGGAGAGTCACC  
TGTGTGTGCAGGGGACACTGTGCTCCTGTGAGAGCTTCATCATGGAGTCAGCCAACGGGTGGCAGCATG  
CATGTTCCAGCTGGAGTTCTTCTGCCCTGACAATCATCTTGTCTGCTCGGTCAACGTTGTTGGAGC  
CTGAGACGGAGGCAGCAGCTGACCAGACAGGCTCGGATGAGGAGGGCCACCCGGTTCATCATGGTGGTGG  
CTTCTGTGTTTCACAGTGTACCTGCCAGCGTGTGGCTAGGCTCTACTTCTCTGGACGGTGGCCAC  
TAGTGCCTGTGACCCCTGTCCACACAGCCCTCCACGTACCCCTGAGCTTACCTACCTGAACAGTATG  
CTGGATCCCCTTGATATTACTTCTCAAGCCCTCGCTCCCAATTCTACCAAGCTCACAATCTGCA  
GCCTGAAGCCCAAACGCCAGGACGCACGAAGACGCGGAGGTGAGAAGAGATGCCAATTTGGAACCTCTG  
CAGTAAGAGCTCCATCGATGGGGCAAATCGTTCCAGAGGCCATCTGACGGGCAGTGGGATCTCCAAGTG  
TGT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online >](#)

<b>ACCN:</b>	NM_175520
<b>Insert Size:</b>	1056 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_175520.4</a> , <a href="#">NP_780729.2</a>
<b>RefSeq Size:</b>	3285 bp
<b>RefSeq ORF:</b>	1056 bp
<b>Locus ID:</b>	243270
<b>UniProt ID:</b>	<a href="#">Q8C131</a>
<b>Cytogenetics:</b>	5 F
<b>Gene Summary:</b>	Acts as a receptor for L-lactate and mediates its anti-lipolytic effect through a G(i)-protein-mediated pathway.[UniProtKB/Swiss-Prot Function]