

Product datasheet for **MC213218**

Eda2r (NM_001161432) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Eda2r (NM_001161432) Mouse Untagged Clone
Tag: Tag Free
Symbol: Eda2r
Synonyms: 9430060M22Rik; TNFRSF27; Xedar
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC213218 representing NM_001161432
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGGATTGTCAAGAGAATGAGTACCGGGACCAATGGGGCGCTGTGTACCTGCCAACAATGTGGCCCTG
GACAGGAGCTCTCGAAGGATTGTGGTTATGGAGAAGGTGGAGATGCACACTGCATAGTCTGCCCTCCCCG
AAAGTATAAAAGCACTTGGGGACATCACAGATGTCAGACATGCATCACTTGTGCTGTCATCAATCGGGTC
CAGAAGGCCAACTGCACAAATACCTCTAATGCTATCTGTGGAGACTGTCTGCCCAGGTTCTACCGAAAGA
CACGCATTGGGGCCCTGCAGGACCAAGAATGCATCCCATGTACAAAGCAGACTCCTTCTCCGAGGTTCA
GTGTACCTTCCAGTTGAGCTTAGTGAAGGTAGATGCACACACTGTTCCCCTAGGGAAGCCACACTTGT
GCACTGGTGGGAAGTCTGCTAGTGGTATTTGCTCTGGCCTTCTTGGACTCTTCTTCTTTATTGCAAGC
AGATCTTCAACAGACATTGTCAATGCAGAGATTCGTTGCAATATGAGGCTGAAAAGACAGTGGAGGAGGA
TTCTCTTTTCCCGTACCACCTGGCCAGGAGACCAGTCTGAATTTCCAGCTAATGAGGGCATCTTGGAG
ATCAAACCACTTAACTCCATCCTGGATGATGACTGCAGCTCCACTCGCGGCTTCCCACCCAGGAGTCTCT
TTACCATGGCCTCCTGTGCCTCAGAGAGCCATTCTCAATGGGTCCATACCCCAATTGAATGTACAGAGCT
GGACCTGCAAAAAGTTTTCCAGCTCTATCCCTCTACTGGACCTGAAACCTTGAGAGAGAACACAGCTGAA
CACTCAGGAGACAGGCTGGAGCTCTATGTACCTTTTGAAGTTCCAGCCTT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001161432
Insert Size: 894 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001161432.1](#), [NP_001154904.1](#)

RefSeq Size: 4798 bp

RefSeq ORF: 894 bp

Locus ID: 245527

UniProt ID: [Q8BX35](#)

Cytogenetics: X C3

Gene Summary: Receptor for EDA isoform A2, but not for EDA isoform A1. Mediates the activation of the NF-kappa-B and JNK pathways. Activation seems to be mediated by binding to TRAF3 and TRAF6 (By similarity).[UniProtKB/Swiss-Prot Function]
 Transcript Variant: This variant (1) represents the longest transcript. All three variants encode the same protein.