

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_139224.1](#), [NP_631970.1](#)

RefSeq Size: 267 bp

RefSeq ORF: 192 bp

Locus ID: 246084

UniProt ID: [Q8R2I3](#)

Cytogenetics: 8 A2

Gene Summary: Has antibacterial activity.[UniProtKB/Swiss-Prot Function]