

## Product datasheet for MC213272

### Olfr1290 (NM\_146262) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Olfr1290 (NM\_146262) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Olfr1290  
**Synonyms:** MOR248-9  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC213272 representing NM\_146262  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGATGAAGGCAATCACACTATGGTGTCAGAATTTATACTCTGGGGACTTTCCAGCTCACAGAGGATTG  
AAGTCTTACTCTTCATGGTGTTCATGCTTTTCAATGCTTTATCTTCTCATTGTTTCTGGAAATATTGTCATTCTGGT  
TTAATCACCACCTGACCCCATCTGCATTCTCCCATGTACTTCTTGTGGCCAACCTGTCCTTTATTGAT  
ATGTGGCTTTCTCAGTCACTACTCCTAAGATGATTACAGACTTTCTCAGAGAAAACAAGACCATTTCTC  
TTGCAGGATGCATGTCCAGGTCTTCTTTGCCATTGTATTGCTGCAGGAGAGATGGTGTCTGTTGGTGTAT  
AATGGCTTATGATCGCTATGTGGCCATCTGCAAACCACTCCACTACTTACCATCATGAATATGAAGAGA  
TGCACTGGGTGGTGTTCACCTTCTGGACAACCTGGATTGTACATGCTCTTAGTCATTTGGTAGTGATTG  
TGGAAATGGCTTTTTGTGGACCTAAGGAAATAGATAGCTTTTTCTGTGACATGCCATTGGTAATCAAGCT  
AGCTTGCATAGATTCCCGTGACTTGGATGTTCTAATGAATGCTGATTGTGGACTTGTGGCTGTAACCTGC  
TTTATTCTGTTGCTTATTTCTACACATATATTCTTATCACCGTATGTCAGAGCTCTAAAGCTGGTGCAT  
CTAAGGCAATGAACACATGCACCTGCCACATTACGGTGGTGTGATCTTTTTGTGCCCTGCATCTTCAT  
CTATGTATGGCCTCTCAATATCACTTGGTTGGATAAATTTCTGTGATTTTACTCTGTTTTTACACCT  
CTCCTAAATCCAGCCATTTATACTCTGAGAAATAAGAGATGAAAAATGCCATGAAAAGATTTATAAGTA  
ACTACCTTAGTCACAAAGGAAATTTA**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_146262  
**Insert Size:** 939 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_146262.3</a> , <a href="#">NP_666374.2</a>
<b>RefSeq Size:</b>	1234 bp
<b>RefSeq ORF:</b>	939 bp
<b>Locus ID:</b>	257662
<b>Cytogenetics:</b>	2 E3
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (1) represents the longest transcript. Variants 1, 2, and 3, encode the same protein. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.