

Product datasheet for MC213279

Olfr1116 (NM_001011734) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Olfr1116 (NM_001011734) Mouse Untagged Clone
Tag: Tag Free
Symbol: Olfr1116
Synonyms: MOR264-8P; MOR264-21P; Olfr1116-ps; Olfr1521-ps1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC213279 representing NM_001011734
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAATTCATTCTTCTGGGTTTTCTAATGTTCCCTATTTGCAGTGGTGCTGTTTATGGTATTTTTGT
 TTATGTATATGACAATTCTGTTGTGCAACAGCATTATAATAGTGTAGCAAAAACAGCCCTGCTCTCCA
 GACTCCTATGTATTTTTCTTAGCAACTTTCTTTTTGGAAATCTGTTATGTAACAGCTACTATTCCA
 AGAATGCTTATGGATTTATATACCCTAAAAGGAAACATTTCTGTGTTGCCTGTGCAACACAATGTATT
 TTGTCCTTACGTTGGGAGGCACAGAGTGCCCTTGTCTAGCAGCTATGGCCTATGATCGTTATGTGGCTAT
 CTGCCACCCTCTCCAATATTCTTTACTTATGAAGAATAAAGTGTGTTTACAGCTGGTAGCTGCCCTCTGG
 ATCAGTGGGATCCAGTAGAAAATGGGCAGACCTACCAGATATTCTCCCTTCACTTTTGTGCTTCTAACA
 GAATAGATCACTTTTTCTGTGACATTCCCCACTTCTCAAGCTTGCTTGTGGTGACACCTTTATGAATAC
 AGTTGCAGTCTATGTTGTTGCAGTGTGTTTGTATGGTTCCATTTCTGCTAATCATTGTCTCTACATC
 AAGATTATCTGCAACATTATGAACTATCTTCAGCCAAAGGGATGGCCAAGGCTTTTTCCACCTGCTCGT
 CTCACCTGATAGTTGTAATCTTGTCTATGGAACAGCAAGCATTACTTACTTACAGCCAAACAAAGTCA
 GTCAGAAGGAATGGGAAGATGCTGTCTTTTTCTACACATTTTGATCCAGCTTTGAATCCTATTATA
 TACAGTCTGAGGAACAAAGATATCATGATGGCACTGAGAAAATTACACAGTAAGTTACTGATATGGTGGG
 AAAACTTAAATAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001011734
Insert Size: 924 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001011734.1</u> , <u>NP_001011734.1</u>
RefSeq Size:	924 bp
RefSeq ORF:	924 bp
Locus ID:	257875
Cytogenetics:	2 D
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]