

## Product datasheet for MC213371

### Olfr967 (NM\_001011826) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Olfr967 (NM\_001011826) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Olfr967  
**Synonyms:** MOR171-30P  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC213371 representing NM\_001011826  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGC**C

ATGATAGCAGGAACTACTCCATGGTGACTGAGTTTATCCTTGCTGGATTAACAAGCACACCAGAAGCTGC  
AGCTGCCTCTCTTCTTCTTCTTCTCGGAATCTATGCAGTAACGATGGTAGGAACTGGGCATGATCAC  
ACTGATTCGCTCAGCTCCACCTGCACACCCCATGTACTTCTTCTCAGCAGTCTGTCCTTCATTGAC  
CTCTGCCATTCAACTGTCATTACCCCCAAAATGCTGGTGAAGTGGTACTGTGAAGAATCATCTCCT  
ACCTGAATGTATGACTCAGCTCTATTTCTTTCTGGTTTTGTTATATCAGAATGTCACATGCTGGCAGC  
TATGGCATATGACCGATACGTTGCCATCTGTAACCCCTTGCTATAACAATGCTATGATGTCCTATCAAGTC  
TGCACCTGGATGATATTTGGGGTGTATAGCATGGGTTTTATTGGTGCCACAGCTCACACAGTGTGCATGC  
TAAGAGTCCATTTCTGTAAGGTTGATGTAATAATCATTACTTCTGTGATCTTTTCCCGTTATTGGAGCT  
TTCTTGCTCCCCTACTTTTATCAATGAGGTAGTTGTTCTATGTTTTAGTGCTTTCAATATCCTTTTCCCA  
ACACTGAGTATCCTGAGCTCTTACATCTTATCATTGCCAGCATCCTCCGATTAATAATCCACTGAAGGCA  
GGTCCAAAGCCTTCAGCACCTGCAGCTCACACATATCAGCTGTTGCTGTGTTCTTTGGGCTGTCGATT  
TATGTACCTACAGCCATCATCAGTCAGCTCCATGGACCAAGGAAAGTCTCCTCTGTATTTTATACCATT  
GTTGTGCCATGCTAAACCCCTTGATCTACAGCCTGAGAAATAAGGATGTCAAGGTCGCACTAACAAAAT  
TTTATGAAAAAGTTTCTCATGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_001011826  
**Insert Size:** 933 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001011826.1</a></u> , <u><a href="#">NP_001011826.1</a></u>
<b>RefSeq Size:</b>	933 bp
<b>RefSeq ORF:</b>	933 bp
<b>Locus ID:</b>	258086
<b>Cytogenetics:</b>	9 A5.1
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]