

## Product datasheet for **MC213454**

### Olfr370 (NM\_146270) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Olfr370 (NM\_146270) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Olfr370  
**Synonyms:** MOR267-16; MOR26716  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >MC213454 representing NM\_146270  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGAGTGTGTAATGACACTGTGGTCAGAGAATTTGTCTTCCTTGGCTTCTCGTCTCTGGCTGAGCTGC  
 AGCTGCTACTCTTTGCCATCTTCTGTCCCTCTATTTGTTCACTCTGAGCACCAACGCTGTAATTGTTTC  
 CACCATTTGTGCTAGACAGAGCCCTTCATACGCCATGACTTCTCCTTTCTGCTCTGCTGTTCTGAG  
 ACCTGCTACACCTTCGTCATTGTGCCAAGATGCTGGTTGACTTGTGGCGGGAAGAAGAGCATCTCTT  
 TCCTCGGCTGCGCATCAAATGTTACCTTCTCCTCGGCTGTTCTCATTCTTCTGCTGGCAGC  
 CATGGGTTATGATCGATACGTGGCCATTTGCCACCCTCTGCGCTACACGGTGTCTATGGGCACAGGGTA  
 TGGTGGGGCTAGTAGCTGCTGCGTGTGTCTGTGGCTTCACTGTGGCACAGGTAATCACATCCCTGGTGT  
 TTCGTCTACCTTCCGCTCTTCCAATCAACTCCACCCTTTTTCTGTGACATCTCCCTGTTCTCCAGTT  
 GGCATCTCACCACCCTCACTCCACTCAGATCACCATCTTCTGCTTTGCGCGTGGTCTTGGTTATCCCT  
 CTCTTCTTGATCCTGGTATCATATTCACATCATTTCTACCATCCTCCAATCCCTTCCACGCTGGGCA  
 GGTACAAAGCGTTTTCCACCTGTGCCTCTCACCTATTGTTGTACATAGTGCATATGGATGTGCCTTTT  
 TATCTACCTAAGGCCAAGTCCAGTACTCTTCCAGCCAGGATGCGCTGATATCAGTATCCTACACTATC  
 CTAACTCCGTTATTCAATCCAGTCATCTACAGCTTAAGAAACAAAGACTTCAAATCAGCTCTTATAGGG  
 TAATAGGAAGAACAGTTACCCTAAGACAACACTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_146270



[View online »](#)

|                               |  |
|-------------------------------|--|
| <b>Insert Size:</b>           | 945 bp   |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).   |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).   |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>  |
| <b>RefSeq:</b>                | <u><a href="#">NM_146270.2</a></u> , <u><a href="#">NP_666382.2</a></u>  |
| <b>RefSeq Size:</b>           | 1096 bp  |
| <b>RefSeq ORF:</b>            | 945 bp   |
| <b>Locus ID:</b>              | 258267   |
| <b>Cytogenetics:</b>          | 8 C2   |
| <b>Gene Summary:</b>          | Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008] |