

Product datasheet for MC213459

Olfr1402 (NM_146275) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Olfr1402 (NM_146275) Mouse Untagged Clone
 Tag: Tag Free
 Symbol: Olfr1402
 Synonyms: MOR217-1
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 Fully Sequenced ORF: >MC213459 representing NM_146275
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGC**C

ATGACTCCATTTGAAATGACCAATCACACACGTGTCACAGAATTCATCTTCTGGGATTTCCAACCACC
 CCAACCTACAGGGTGTCTTCTTCTGGCCTTCTGGCCATTTACCTGACAACCTCTCTGGGAACACACT
 CATGATAGTGGCCACCAGGGTCAGCCCTGCACTCCACACTCCAATGTAATTTCTCAGCAACTTGAGC
 TTCTGGACATCTGCTACACATCTACCTCCATTCCAGTCATGCTAGTGAATTTCTCAGGAGAAGAAGA
 CCATCTCCTTTGAGGGCTGCCTCTCTCAGATCTTCTTCTTTGTGTCATGTGCTGGAAGTGAATGTGCTT
 GTTGGCTGCCATGGCATAATGACCGCTATGTGGCCGTTTGTACCCTCTTCGGTATCCAGTTCTCATGAGC
 ACAAGGGTCTGTATCAGCTTGGTACTGGATCCTGGCTCTGTGGCTTGGTGAATTTCTGTGACTCACACAG
 TATTGACAAGCACACTTACTCTGTGTGGCCCAACCAATCAGCCACTTCCTCTGTGATGTCCCATTGCT
 CCTGAAGCTCTCCTGCTCAGACACCTCTGTCAATGAATCTGTCTCCATGTGTCCAGTGCCACTGTTGGC
 CTGAGTCCCTGCCTGTTCACTGCAGGCTCCTATATACTGATCATCTCTGCCATCCTGAGGATCCCTCTA
 CTCAGGGCAAGAGAAAGGCCTTCTCCACCTGTGCCTCACACCTCACTGTGGTGGTGGTTTTCTTTGGAAC
 AGCCAACCTCAACTATGTGAGCCCAAGGAAGGCTACTCCCTGGACATGGGCATATTAGTGTCTGTGCTC
 TACTGTGTTGTGACCCCTGTTAATCCATCATCTACAGCCTGAGGAACAAGGAGTCAAAGGTGCGC
 TGAGAAAAGTACTGAGTGCCTTCTTAAACAATAGCAGTGGCTAGAGAG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
 ACCN: NM_146275
 Insert Size: 963 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_146275.1</u> , <u>NP_666387.1</u>
RefSeq Size:	963 bp
RefSeq ORF:	963 bp
Locus ID:	258272
Cytogenetics:	3 F2.1
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]