

## Product datasheet for MC213560

### Olfr1284 (NM\_146381) Mouse Untagged Clone

#### Product data:

|                           |                                                                                 |
|---------------------------|---------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                             |
| Product Name:             | Olfr1284 (NM_146381) Mouse Untagged Clone                                       |
| Tag:                      | Tag Free                                                                        |
| Symbol:                   | Olfr1284                                                                        |
| Synonyms:                 | MOR245-13                                                                       |
| Mammalian Cell Selection: | Neomycin                                                                        |
| Vector:                   | pCMV6-Entry (PS100001)                                                          |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                            |
| Fully Sequenced ORF:      | >MC213560 representing NM_146381<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGATCAAATAAATGAACTGTAGCTTTTGAATTTGTGTTGTTGGGACTCTCCAGTTCTTGAAAAATA  
CAATTTTCTCATGTCCACCTTCTTTGCTCTATGTAAGCATCATCGTGGGAAACCTTTTCATTGTCTT  
TTTGGAATTAATGACTCCCGTTTACAGTCTCCTATGTATTTCTTTGGCCAATCTGTCCTCATTGAT  
GTTGGACTTTCTCTACCACAGTGCCCAAGATGATCTCAGATCTTCTGAAAGAGCACAAAGTAATATCTT  
TCCACAGTTGCATGACTCAGATCTGCTCCATCCACATTATGGGAGGAGTGGAGATGGTGTGCTCATAGC  
CATGGCGTTTGACAGGTACACGGCCATCTGTAAGCCTCTGCACTACATGAGCATCATGAGCCCTAGAATA  
TGCAATTTCAATTTGTAATTGCTGGCTGGTCACTGGAGTGGTCCATGCCATGTACAGTTCTCTTTTGTG  
TGAACCTGCCTTTTGTGGTCTTAACAAAGTGGACAGCTTTTACTGCGACTTCCCTCGAATCATACAAC  
GGCATGCACGGATAGAGACACATTTGAGTTTGTGTTGCTGCCAACAGTGGCTTCATGACCTTGGGGACC  
TTCTTCTTGCTTCTCCTCTCCTATGTCTTATTTTGGTCACAGTCTGGCAAAGGTCTTCAGGGGACTTGT  
CAAAGGCGCTTGCTACTCTGTCACTCACATCACTGTGGTTGTCTGTTTTTACTCCATGCATGTTTCT  
CTATGTGTGGCCTTTCCACATCATCGACTGACAAATACCTGTTCAATGTTGACTTTGCTGTACCCCT  
GCTCTTAATCCTGTCACTATACTTTAAGGAACAAAGATATGAAGGAAGCCATAAAAAAGCTGAGCAAGC  
AGAGATGTTACATCAGAATTTTC**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

|                    |           |
|--------------------|-----------|
| Restriction Sites: | SgfI-MluI |
| ACCN:              | NM_146381 |


[View online »](#)

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Insert Size:</b>           | 936 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>                                                                                                                                                                                                        |
| <b>Note:</b>                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>RefSeq:</b>                | <a href="#">NM_146381.1</a> , <a href="#">NP_666493.1</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>RefSeq Size:</b>           | 936 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>RefSeq ORF:</b>            | 936 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Locus ID:</b>              | 258379                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Cytogenetics:</b>          | 2 E3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Gene Summary:</b>          | Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008] |