

## Product datasheet for MC213591

### Olfr291 (NM\_146415) Mouse Untagged Clone

#### Product data:

Product Type: Expression Plasmids  
Product Name: Olfr291 (NM\_146415) Mouse Untagged Clone  
Tag: Tag Free  
Symbol: Olfr291  
Synonyms: MOR254-2  
Vector: pCMV6-Entry (PS100001)  
E. coli Selection: Kanamycin (25 ug/mL)  
Cell Selection: Neomycin  
Fully Sequenced ORF: >MC213591 representing NM\_146415  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGAGATGGACGTCTACAATCTTACCACAGTGACTCAGTTCATCCTCATAGGGCTCTGACCTCCCTG  
AAGTGCCTATCCACTCTTTGTGGCCTTTGTGCATCATCTATCAGATCACTTTACTGGGAAATGGGCTCAT  
CCTCTTGGCCATTGTGACTGAGAAAAAGCTTCAAACCTCCCATGTATTACTTGTGGCAAATCTGTCCTTA  
CTGGACATATTCTGCCATCAGCTACTGTCCCAAGATGCTCAAGAATCTCTTGACTGAGGATCACAGCA  
TTTCCTTTGTTGGGTGTGCTTTACAGCTCTATTTCTGGTGGCCCTAGCTGGGACTGAAGCTTCTTGCT  
GGCTGTGATGGCTTATGACCGGTATGTGGCCATATGCTTTCTCTACGTTACTCTCTCATTATGACCAAG  
GTTTCGCTGTGCAGCTGTTGTTGGGACTTGGGCAGCTGGGTTTCTGAACTCTTTGTCCACACAATGT  
CCACCTTAGCCTTTCTTTCTGCAAGTCTAATAGAGTTAACCAGTATTATTGTGATATCCACCTGTGGT  
GGCCCTGTGATGCTCATCCACCTATATGGCAGAAGTGCTTCTTTAGTGATAGCAAGTATCTTTGGGGTT  
GGTGCTTTCTGATCACTTTGATCTCTACATATACATTGTCTCCACCATCCTAAAGATCCGGTCAGCAG  
AAGGAAAGCGCAAAGCTTCTCCACATGTGCTTCCCATCTTCTGTAGTCTTCTGTTCTATGGTACAAC  
TATATTTACCTATATGCGCCCAACCTCCAGTCAACACTCTCTGGAAGAGACAGACTCATCTCTATGTTG  
TATGGGGTCACTTACTCCCATGTTAAACCCATTATCTACAGCCTGAGAAACACAGAAGTCAAAGGAGCTC  
TAAGAAAAGTTTTACATCTTCGGATATGTTACAGACAGCATGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI  
ACCN: NM\_146415  
Insert Size: 954 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_146415.2</a> , <a href="#">NP_666527.1</a>
<b>RefSeq Size:</b>	1159 bp
<b>RefSeq ORF:</b>	954 bp
<b>Locus ID:</b>	258410
<b>Cytogenetics:</b>	7 D3
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]