

Product datasheet for MC213598

Olfr470 (NM_146425) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Olfr470 (NM_146425) Mouse Untagged Clone
 Tag: Tag Free
 Symbol: Olfr470
 Synonyms: MOR204-22
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 Fully Sequenced ORF: >MC213598 representing NM_146425
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGCTTTCCTGGAGGATGGGAACCACTGCAGTGACAGAGTTCGTTTTATTTGGATTAAACAGATGATC
 CAGTTCCTAGAGTCATCCTCTTCATCATCTTCTGTGCATCTACCTGGTGAATGTGTCTGGGAACCTCAG
 CACCATCCTTCTCATCAGAGTCTTCCAGCTCCATCACCCATGTACTTTTTCTCAGCCACTGGCT
 TCTGTTGACGTAGGCTACTCATCTACTGTACACCCAAAATGCTTGCCAACTTCCCTGCTGGAGAGAAGCA
 CAATCTCTACCTTGGGTGTACCATCCAGCTTTTCTCAGGTGCTTTTGTGGGACACTTGAGTGCTTCT
 TCTGGCCACCATGGCTTATGATCGTTTTATAGCAATCTGCAACCCATTGCTTTATCCACCAAAATGTCC
 ACACAAGTCTGTATCCAGTTGCTTGTAGGATCTTATATTGGTGGGTTTCTTAATGCTTCTCCTTTCTCC
 TTTCTTCTTCTCTTCTTCTGTGGACCAATAGAGTCAATCACTATTCCTGTGATTTGACTCCTTT
 GATAGAGCTCTCCTGTTCCAGGTAGCAATGTCCCATAGTTCCTGCCTCATTTTGTCTGCCTTTGTCTT
 ATAGTACAGTGTCTGTATAGCCATCTCTACACCTACATCCTCATCACCATCCTGAAGATGCGCTCCA
 CTGAGGGCCGACAGAAGGCCTTCTCCACCTGCCTTCCACCTCACTGCAGTGACTCTCTATTATGGAAC
 TGTACATTCAATTTATGTGATGCCAAAATCCAGCTACTCCAGGACCAGAACAAAGGTGGTGTCTGTGTT
 TACACAGTGGTAATCCCAATGTTGAACCCATCATCTACAGCCTCAGGAACAATGAGATTAAAGGTGCTC
 TGAAGAGACAGCTTGCTAGGAAAATATTCTCT**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
 ACCN: NM_146425
 Insert Size: 945 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_146425.1 , NP_666636.1
RefSeq Size:	945 bp
RefSeq ORF:	945 bp
Locus ID:	258417
UniProt ID:	Q8VF65
Cytogenetics:	7 E3
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]