

Product datasheet for MC213682

Olfr99 (NM_146515) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Olfr99

Synonyms: MOR156-1

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC213682 representing NM_146515
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAACTGCAGTCAGGCTCCTACATTAATCCTCTTGGGACTGTCTAGTGATGCAGAGAAATGGCAGCCCC
TCTTTAGCATCTTCCTTGTTATCTACTTGCTGGGTCTCTTAGGGAACCTGCTACTTTTGTTAGCTATTGG
CACTGATGTCCACCTCCACACCCCATGTATTTCTTCTCAGTCAGCTCTCATTTGTGGATCTTTGCTTC
ATCACCACCACAGCCCCATAAATGCTGGAGGCTCTGTGGACTGGAGATGGATCAATCTCATTCTCTGGAT
GCCTAACTCAGTTGACTTCTTTGCTGTTTTGCTGACATGGATAACCTGCTTCTAGCAGTCATGGCTAT
TGACAGGTATGCCGCTATCTGTCATCCACTGCGCTACTCAGCTCTCATGACTCCTTTTAGATGTGGGGTT
CTAGTCAGTGGGTATGGGGAGTAATAATTGTGTCTCTGACCCATACTTTGTTGCTCTCTAAGTTAT
ATTTCCATACCAATCAAGAGATTCCTCATTTTTCTGTGAATTTGGGCTCTCTATTGCTTTCTGCTC
AGATACTCACCTCAATAAAATCCTGGTGATAATTTGGTTGGAATTTAGGAATCAGTGCGGTTCTCTGC
ATTGTAAGTTCTTATGTTGTATTTTTATGCTGTGGCTAAGGTTCCATCAGCACAGGGGAAAAGGAAAG
CATTGTCCACGTGCAGTCCCACCTCTCTGTAGTCCTCTTTCTACAGCACAGTCTTTGCTACCTACCT
GAAGCCCCCATCTAGTTCTCGTTCTCTGAGGAGGTGGTAGCTGCTGTATGATTCTCTCGTAACCTCC
ACTCTGAATCCCTTCATTTATAGCTGAGAAATAAGGATGTTAAAGTTTCATTGAGGAGAATCCTAAACA
TGGAGTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_146515



Insert Size:	918 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_146515.2</u> , <u>NP_666726.1</u>
RefSeq Size:	918 bp
RefSeq ORF:	918 bp
Locus ID:	258508
Cytogenetics:	17 B1
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]