

## Product datasheet for MC213691

### Olfr860 (NM\_146528) Mouse Untagged Clone

#### Product data:

Product Type: Expression Plasmids  
 Product Name: Olfr860 (NM\_146528) Mouse Untagged Clone  
 Tag: Tag Free  
 Symbol: Olfr860  
 Synonyms: MOR146-2  
 Vector: pCMV6-Entry (PS100001)  
 E. coli Selection: Kanamycin (25 ug/mL)  
 Cell Selection: Neomycin  
 Fully Sequenced ORF: >MC213691 representing NM\_146528  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGAAGTATATAATTTAACAAGTAACTTGAATTCTCTCTCTGGACTCTCAGAGGACCCCTGAAGTAC  
 AACCTGTCTTATTTGCACTATTCCTGCTCATATACCTGCTCACTGTGCTTGGAAATGTGCTCATCATCTT  
 GGCCATCAGCTGTGATTCCACCTCCACAGCCCAATGTACTTCTCTCTACAACCTTTTCATTGTCTGAC  
 ATGGGCTTCAGCAGCACCACCATCCCCAAAATGCTGATAAACTTGCATGCCATAACAGAACCATAACTT  
 ATGCTGAATGCCTAACCCAGGTATCTTTCTTTTTCTTTTTGGGTGTATGGACAGCGTTCTACTGGCAGT  
 GATGGCCACGACCGATGGGTGGCCATTTGTCACCCTCTACACTACCAAGTCATTCTAAATCCTCGTCTG  
 TGTAGATATTTGGTTGTAATGTCATTTTGTATCAGTCTCATTGATTACAGGTGCACTGCTTTATGGTGT  
 CACAATAAAAATTTGTAATAATAAAAAATCCCTCATTCTCTGTGATGTTCCAGAGCTTCTAAAAT  
 TGCTGTCTGACACCTATAAATAGCATAGTTATTTCTTGTGATCATTGTTGGATTCCCTCCCT  
 GCCTCAGGAATCTTTACTCTACTATAAGATTATTTCTCCATTGTTAGAGTTCCTTCACTCTGGGGGA  
 AATGTAAGCGTTCTCTACTTGTGGATCTCATCTGTCAGTTGTTGCTTATTTTATGGAACAGGCTTGG  
 TGTATACCTTAGCTCATCTATTTCTAGTTCTTCCAAGGAAAGTGGTGACTTCAGTAATGTATACCATG  
 GTGGTTCCATGATAAACCCCTCATCTACAGCTTGGGAACAAGGATATCAAAAAGCCCTACAGAAAA  
 TTTTCAGTCAAATAATAATGTTACCTACCTATATAATCCCT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI  
 ACCN: NM\_146528  
 Insert Size: 954 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_146528.2</a> , <a href="#">NP_666739.1</a>
<b>RefSeq Size:</b>	1219 bp
<b>RefSeq ORF:</b>	954 bp
<b>Locus ID:</b>	258521
<b>Cytogenetics:</b>	9 A3
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]