

## Product datasheet for MC213694

### Olfr1168 (NM\_146531) Mouse Untagged Clone

#### Product data:

Product Type: Expression Plasmids  
 Product Name: Olfr1168 (NM\_146531) Mouse Untagged Clone  
 Tag: Tag Free  
 Symbol: Olfr1168  
 Synonyms: MOR174-13  
 Vector: pCMV6-Entry (PS100001)  
 E. coli Selection: Kanamycin (25 ug/mL)  
 Cell Selection: Neomycin  
 Fully Sequenced ORF: >MC213694 representing NM\_146531  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGATAAACCAGATTCTGTAACCACATTTATTCTGGTGGGTTTCTCAGAATATCCACACCTACAGCCAC  
 CCCTCTTCCTTATGATAATGACCATTTACACAGTCACTCTGGTAGGAAATGTAGGCATAATTTTGGTCAG  
 AAGGATCAACCCTAAGCTTCACACACCCATGTACTTTTTCTGAGCCATCTTTCATTTTTGGATATATGT  
 TATTCTAGTGTATTTACACCAAATGTTAGAGATCTTGATTGTGGAACATAGAACTATTTCTTTAAATG  
 GTTGATGACACAGTTTTCTTTGGCTGTGCATGTGTGATTACAGAGATGTTTCATGTTAGCAGTGATGGC  
 CTATGACAGGTTTGTGGCTGTTGCAACCCCTGCTCTATACTGTTGCTATGTCTCATCAATTGTGTGCT  
 CTCCTAGTAGCTGGAAGTTATATGTGGGGAGGTCTGTGTGCTGTGATAATCACATACACTCTGTAGAGC  
 TCTCTTATTGTGAACCTGGCATCATTGATCACTTTGGCTGTGAATATTCAGCCATTGTCTGTATCCTG  
 TTCTGACCCCTCCTTCAGCCAGATGGTATGCTTAGTCATTTCCATATTAAGTGAGGGTAGTAGCCTTCTT  
 ATCACTATGGCCTCCTATGTCTTCATCGTTGTCAATCATCAAGATGCCTTCTAAAGATGGACTCAGAA  
 AAGCCTTCTCCACTGTACTTCACACCTGACTGCCATCAGCATCTCCATGGGATTATTCTCCTTCTCTA  
 CTGTATACCCAATGCAAAGAGCTCAAAGCTGCTGGTCAAAGTAGCAACTGTACTCTATACAGTCCTTATA  
 CCCATGCTGAACCCCTCATCTACAGCCTAAGGAACAAAGATGTAAAAGAGACAGTCAAGAGGCTCATCA  
 GCTCAAAATTGCACTCTCAGACAATATAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI  
 ACCN: NM\_146531  
 Insert Size: 939 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_146531.2</a></u> , <u><a href="#">NP_666742.2</a></u>
<b>RefSeq Size:</b>	939 bp
<b>RefSeq ORF:</b>	939 bp
<b>Locus ID:</b>	258524
<b>Cytogenetics:</b>	2 D
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]