

## Product datasheet for MC213739

### Olfr1020 (NM\_146580) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr1020 (NM_146580) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Olfr1020
Synonyms:	MOR201-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC213739 representing NM_146580 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGTCAGGAGTGGAAAAGGAATCCAAATAAGAATGCCACAGAAGTGACAGAATTTATCCTTTTGGGTC  
TCTCAGACAATCCAGACTTGCAAGGTGTCCTCTTTGCATTATTCCTGATCATCTACAGGATGACTGGT  
GGGTAATTTGGGCATGATGGCCCTGATTAAGATTGACCGCAGTCTGCACACACCTATGTACTTCTTTCTC  
AGCAGCCTTTCTTTGTTGATGCCTCTTATTCTTCTGTCAACCCCAAGATGCTAGTGAACCTCATGG  
CTGAGGATAAGAGCATTCTTTCAATGGATGTGCTACCCAGTTCTTTTCTCGGTTCTTCTTGGGAAC  
TGAATGCTTCCTGCTAGCCATGATGGCCTATGACCGCTATGCAGCCATCTGGAATCCCCTGCTATACCCA  
GTTCTCATGTCTGGGAGAATTTGTTTCATGTTGGTGTCTACTTCATTCTTAGCAGGCTTTGGCAATGCAG  
CCATCCACACAGGGATGACTTTTCACTGTCTTTTGTGGCTCTAATAAGATCAATCATTCTATTGTGA  
CACCCACCCCTGCTCAAACCTCTTGTCTGACACTCACATTAATGGCATTGTGATCATGGCCTTCTCC  
AGTTTCAATGTCATCAGTTGTGTTCTGATTGTTTTGATATCTTACTTGTGCATCCTCATTGCCATACTGA  
AGATGCCTTCAGCAGAGGAAGACACAAAGCCTTCTCCACCTGTGCCTCCACCTCATGGCTGTTACCAT  
CTTTTTTGAACAATTCTTTTCATGTAATGCGCCCTACCTCTAGCTACTCTATGGAACAGGACAAAGTT  
GTCTCTGTGTTTATACAGTAGTCATTCCAATGCTTAATCCTCTCATCTATAGTTTAAAAATAAAGATG  
TTAAGAAAGCAGTGAAGAAAATCTTACATAATTATGTGGTT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_146580


[View online »](#)

<b>Insert Size:</b>	954 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>Note:</b>	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
<b>RefSeq:</b>	<a href="#">NM_146580.2</a> , <a href="#">NP_666791.1</a>
<b>RefSeq Size:</b>	1079 bp
<b>RefSeq ORF:</b>	954 bp
<b>Locus ID:</b>	258573
<b>UniProt ID:</b>	<a href="#">Q8VFK7</a>
<b>Cytogenetics:</b>	2 D
<b>Gene Summary:</b>	<p>Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]</p>