

Product datasheet for MC213935

Olfr1248 (NM_146791) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Olfr1248 (NM_146791) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Olfr1248
Synonyms:	MOR23110
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC213935 representing NM_146791 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGACAGAATAACAATGTTACAGAGTTTATCCTCCTGGGTCTCACTCAGGATCCTGCTGGGCAAAAAG
 TATTGTTTGTATGTTTTACTCATCTACATTGTGAAATAGTGGGCAACCTGCTCATTGTGGGACAGT
 GATTGCCAGCCCCTCCTGGGCTCCCCAATGTACTTTTCTTGCCCTTTCTTTCCTCATGGATGCTGTT
 TATTCTACTGCCATCTTGCCCAAGTTGCTTACAGATTTACTTTGTGATAAGAAGACCATCTCCTTCACAG
 CTTGTCTGGTTCAGCTTTTGTGGAGCACTTATTTGGTGGTCTGAGGTCTTCATTTTGGTGGTATGGC
 CTATGATCACTATGTGGCCATCTGTAAGCCACTGCACTATCTGACCATAATGAAACGGCAGGTTTGCATT
 CTCTTGGTGGTGGTGTCTGGGCTGGAGGATTTGCACATGCTCTGCTTCAAGTTATCTCTGTGTATTTAC
 TTCCTTTTGTGGACCCAATGTCATTGACCATTTTGTCTGTGACATATACCCATTGTTAGGACTTGCATG
 TACTGACACTTACTTCTTGGACTCACTGTAGTTAGCAGTTATGGAGCAATGTCTATAGTGGCCTTTATC
 CTCCTTCTGTCTCCTATGGAATCATTCTAACTCTCTTAAGACTGATAGTCAGGAAGGGAGGCGCAAAAG
 CCCTGTCCACCTGCAGCTCCACATCATGGTGGTGTCTCTTTTGTTCCTGTATTTTCATGTATGT
 TAGACCTATCTCCAACCTTCCAATTGATAAATATACTGTCTTTTATACAATTTTACTCCCATGTTG
 AACCCCTTAATATACCTTGAGAAACATGGAGATTAATAAATTGTATGGCAAAGCTGTGGAGTAAATGT
 TCACTAAAGCATAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_146791


[View online »](#)

Insert Size:	924 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_146791.2</u> , <u>NP_667002.2</u>
RefSeq Size:	1078 bp
RefSeq ORF:	924 bp
Locus ID:	258787
Cytogenetics:	2 E1
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]