

Product datasheet for MC214005

Olfr147 (NM_146869) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Olfr147 (NM_146869) Mouse Untagged Clone
Tag: Tag Free
Symbol: Olfr147
Synonyms: M3; MOR164-1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214005 representing NM_146869
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGC**C

ATGATAAGCATGCTGGCTGGAAATGGCTCTTCTGTGACAGAATTTGTCCTTGCTGGTTTGACAGATCGTC
 CAGAGCTCCAGCTGCCTCTCTTTACCTGTTTCTAATAATCTACATAATCACAGTGGTGGAAACTTGGG
 CTTGATCATCCTGATTGGCCTCAATCCTCACCTGCACACCCCATGTAATTTCTCTTCAACCTCTCC
 TTCATTGATCTCTGTTACTCTTCTGTCTTCAGCCCCAAAATGCTGATTAACCTTTGTCTCTGAGAAGAATT
 CCATCTCCTATGCGGGTGCATGACTCAACTGTTTCTCTTTCTCTTTTTGTCATCTCTGAATGCTACAT
 GTTGACCTCAATGGCCTATGATCGCTATGTGGCCATCTGTAATCCACTATTGTATAAGGTCACCATGTCC
 CCTCAGATCTGTTCTGTGATATCTTTGCTGCGTATGGGATGGGATTTGCTGGATCCTCTGCCACACAG
 GCTGTATGCTCAGACTGACCTTCTGCAATGTCAATGTCATCAACCATTACTTATGTGACATTCTTCTCT
 CCTTCACTTTCTGCACCAGTACCTATGTCAATGAGGTTGTGGTTCTCATAGTTGTGGGATTAACATC
 ACAGTTCCAAGCTTACCATCTCATTTCTATGTTTTCTCCTTGCCAACATTCTAAACATCAAATCCA
 CACAAGGAAGAGCAAAGCCTTACGACCTGTAGCTCTCACATCATGGCAATTTCTTTTCTTTGGATC
 AGCTGCATTTATGTATCTTAAATATTCTTCTGGATCTATGGAACAAGGAAAGATATCTTCAGTTTTCTAC
 ACTAATGTTGGTCCCATGCTCAACCCTCTGATTTACAGTTTGAGGAATAAGGATGTCAAGGTGGCATTAA
 GGAAATCATTGATTAAGTTCAGAGAAAAGACAGATTTTAA**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_146869
Insert Size: 954 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_146869.2</u> , <u>NP_667080.1</u>
RefSeq Size:	1175 bp
RefSeq ORF:	954 bp
Locus ID:	258869
UniProt ID:	<u>Q60886</u>
Cytogenetics:	9 A5.1
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]