

Product datasheet for MC214052

Olfr376 (NM_146922) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Olfr376 (NM_146922) Mouse Untagged Clone

Tag: Tag Free

Symbol: Olfr376

Synonyms: MOR135-12

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC214052 representing NM_146922
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCAGGGCACCACAGAAAGGAATCAAACCTGCCATCTCCAGTTCCTCCTCTGGGCCTGCCATCCCCA
CAGAGCATCAACACCTGTTCTATGCCCTGTTCTGGCCATGTACCTCACCACCGTCCTGGGGAACCTCAT
CATCATCATCCTTATTCATCTGGACTCCCATCTCCACACCCATGTATTCTTTCCTCAGCAATTTGTCC
TTCTCTGATCTCTGCTTTTCATCTGTACAATGCCAAATTGCTGCAGAACATGCAGAGCCAGGATCCAT
CAATTCCTATGCAGGCTGCCTGGCACAATGTACTTTTCTGTTTTTGCAGACCTTGAAAGCTTCCT
CCTTGTGGCCATGGCCTATGACCGCTATGTGGCTATCTGCTTCCCCTTCATTACATGAGCATCATGAGC
CCCAGGCTCTGTGTGAGTCTGGTGTCTGCTCCTGGGTGCTGACCACCTTCCATGCCATGCTGCATACCC
TGCTCATGGCCAGATTGTCACTCTGTGAGGACAATGTGATCCCCCACTTTTTCTGTGACATGTCTGCTCT
GCTGAAGCTGTCTGCTCTGACACCTATGTTAATGAATTGGTGATATTTGTATGGGAAGCCTGATCCTT
GTCATTCCATTTGTGCTCATCCTTGTGCTCTATGCACGAATTGTGCTCCTCATTCTCAAGTCCCCTGCTG
CTCGAGGCATCCGTAAAGCCTTCTCCACCTGTGGGTCCACCTGTCTGTGGTGTCACTGTTCTATGGAAC
AGTCATTGGTCTGTACTTATGTCCATCAGCTGATAACTCTACTGTGAAGGAACTGTATGGCCATGATG
TACACAGTGGTGACTCCCATGCTGAACCCCTTCATCTACAGCCTGAGGAACAGAGACATGAAGGGGGCCC
TAATTAGTGTTCTTTGAAGAAGAAAATTCTTTTCTGTCTA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_146922



[View online »](#)

Insert Size:	954 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_146922.2 , NP_667133.1
RefSeq Size:	1079 bp
RefSeq ORF:	954 bp
Locus ID:	258924
Cytogenetics:	11 B4
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (2) differs in the 5' UTR and uses a downstream in-frame start codon, compared to variant 1. The encoded isoform (2) has a shorter N-terminus than isoform 1. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.