

Product datasheet for **MC214053**

Olfr376 (NM_001172686) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Olfr376 (NM_001172686) Mouse Untagged Clone
Tag: Tag Free
Symbol: Olfr376
Synonyms: MOR135-12
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214053 representing NM_001172686
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGCCAGATGCAGGGCACCACAGAAAGGAATCAAACGCCATCTCCCAGTTCCTCCTCTGGGCCTGC
 CCATCCCCACAGAGCATCAACACCTGTTCTATGCCCTGTTCTGGCCATGTACCTACCACCGTCTGGG
 GAACCTCATCATCATCCTTATTCATCTGGACTCCCATCTCCACACCCATGTATTCTTTCCTCAGC
 AATTTGTCTTCTGATCTCTGCTTTTTCATCTGTGACAAATGCCCAAATTGCTGCAGAACATGCAGAGCC
 AGGATCCATCAATTCCCTATGCAGGCTGCCTGGCACAAATGTAATTTTTCTGTTTTTGCAGACCTTGA
 AAGCTTCTCCTTGTGGCCATGGCCTATGACCGCTATGTGGTATCTGCTTCCCCCTTATTACATGAGC
 ATCATGAGCCCCAGGCTCTGTGTGAGTCTGGTGTCTGTCTGGGTGCTGACCACCTCCATGCCATGC
 TGCATACCCTGCTCATGGCCAGATTGTCATTCTGTGAGGACAATGTGATCCCCACTTTTTCTGTGACAT
 GTCTGCTCTGTAAGCTGTCTGCTCTGACACCTATGTTAATGAATTGGTGATTTTGTGATGGGAAGC
 CTGATCCTTGTGATTCCATTTGTGCTCATCCTTGTGCTTATGCACGAATTGTCTCTCCATTCTCAAGG
 TCCCGTCTGCTGAGGCATCCGTAAGCCTTCTCCACCTGTGGTCCCACCTGTCTGTGGTGTCACTGTT
 CTATGGAACAGTCATTGGTCTGTAATGTCATCAGCTGATAACTCTACTGTGAAGGAAACTGTCATG
 GCCATGATGTACACAGTGGTACTCCCATGCTGAACCCCTTCTACTACAGCTGAGGAACAGAGACATGA
 AGGGGGCCCTAATTAGTGTCTTTGCAAGAAGAAAATTCTTTTCTGTCTA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001172686
Insert Size: 963 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001172686.1</u> , <u>NP_001166157.1</u>
RefSeq Size:	1210 bp
RefSeq ORF:	963 bp
Locus ID:	258924
Cytogenetics:	11 B4
Gene Summary:	<p>Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and it encodes the longer isoform (1). The coding region contains an intron, which is atypical for olfactory receptor genes, so translation initiation may actually occur at a downstream in-frame start codon.</p> <p>Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>