

Product datasheet for MC214079

Olfr341 (NM_146950) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Olfr341 (NM_146950) Mouse Untagged Clone

Tag: Tag Free

Symbol: Olfr341

Synonyms: MOR136-2

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC214079 representing NM_146950
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGAGGAGAGATAATGAGAGCACTGTGTCTGAGTTCATCCTCCTGGGGCTCCCCATTAGCCAGAGGATC
 AAGGCCTGTATTCTGCTCTGTTCTGGCCATGTACCTGACAACTGTGCTGGGGAACCTGCTCATCATCCT
 GCTCATCAGGCTGGACTCTCACCTCCACACCCCATGTATTTCTCCTCAGCCACTTGGCCTTCACGGAC
 GTCTCTTTCTCATCAGTCACAGCTCCAAAGATGCTCATGAATATGCTGACACACAGTCAGTCCATCTCAT
 ATGCTGGTTGTGTTTCTCAGGTATACTTTTTTCAACTTTTACTGATCTTGACAGTTTTCTTCTGACCTC
 CATGGCCTATGACAGGTATGTGGCCATCTGCCACCCTCTGCACTATACCACCATCATGAGTCAGAACCTC
 TGTGTCTTTTAGTAGTAATGTCTGGGTTTTATCCTCTGCCAATGCCCTTGTGCACACCCTTCTTTGG
 CTCGCCTATCTATTTTAGAAATAATACCATTCCCCACTACTTCTGTGAACCCTCTGCCTTATTGAGCCT
 GTCGAGCTCAGACACCACCATCAATGAGATGGTCATCCTTCTTTAGGTACACTAGTCATTACCCTGCCA
 TTCATATGCATCCTGGTCTTTATGGCCGAATTGGAGTCACCATCTGAGAACTCCGTCCATCAAGGGAA
 TCTGCAAAGCCTTGTCACGTGTGGTTCTCACCTCTCTGTTGTTGTCTGTAATGAGCTATTATTGG
 GCTTTACCTAGTTCCTCATCTAATAACACCAATGACAAGGATGTCATTGTGGCTGTGATATATAGTCTG
 GTCACACCATGGTGAATCCCTTTATCTATAGTCTGAGGAATCGGGATATAAAAGGAGCATTGAGAAATA
 TCCTCAACAGGAGACTGTGTCCACAGTGG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_146950



[View online »](#)

Insert Size:	942 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_146950.1 , NP_667161.1
RefSeq Size:	942 bp
RefSeq ORF:	942 bp
Locus ID:	258952
Cytogenetics:	2 B
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]