

## Product datasheet for MC214148

### Olfr380 (NM\_147025) Mouse Untagged Clone

#### Product data:

Product Type: Expression Plasmids  
 Product Name: Olfr380 (NM\_147025) Mouse Untagged Clone  
 Tag: Tag Free  
 Symbol: Olfr380  
 Synonyms: MOR135-1  
 Vector: pCMV6-Entry (PS100001)  
 E. coli Selection: Kanamycin (25 ug/mL)  
 Cell Selection: Neomycin  
 Fully Sequenced ORF: >MC214148 representing NM\_147025  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGATAAAGAACAACCAAACTCTCTTCAACCAGTTCTCTCCTGGCCCTGCCGATCCCTGCAGAGCACC  
 AGCAACTGTTCTTTGCCCTGTTCTGTCCATGTACCTCACCACCATCCTGGGAACCTCATCATCATCAT  
 CCTCATTGACTGGACTCCCATCTCCACACCCCATGTACTTGTCTCAGCAACTGTCTCTCTGAT  
 CTCTGCTTTTCTCTGTGACAATGCCAAAGTTGCTGCAGAACATGCAGAGCCAGGACACATCCATCACCT  
 ATGCAGGCTGCCTGACACAAATGTACTTTTTTGTCTTTTTGGAGGCCTGGAAATCTTCTCCTTGTGGT  
 CATGGCCTATGACCGCTATGTGGCCATCTGCCTTCCCCTTCACTACACCAGCATCATGAGTCTCAAGCTC  
 TGTGTGTGTCTGGTGTCTGTCTGGGTGATTAGTATACTGAATTCCATGTTACACACTCTACTCTTGG  
 CTAGGTTGTCTTCTGTGAGGACAACATGATCCGCCACTTTTTCTGTGACATGTCTGCCCTGCTCAAGTT  
 GGCTGCTCTGATATTTATATAAATGAAGTTATGATATTTATCTTGGGAGGACCCCTTATGGTCATTCCA  
 TTCTTACTCATTGTTATGCTCTATGTACAAATTATCTTCTCTATTCTAAAGGCTTCATCTACTCGGGCCA  
 TCTACAAGTCTTCTCTACCTGTGGCTCCACCTGACTGTGGTCACTATTCTATGGGACAATTATTGG  
 TCTCTACTTATGCCATCAGCTAATAACTTTACAGCAAAGGAGGCTTCTATTGCCATAATGTACACAGTG  
 TGACTCCCATGCTAAACCCTTTTCATCTATAGCCTGAGGAACAGAGATATAAAGGAGGCCCTTATAAACG  
 TTCTTATCAAGAAAATACCTTTATAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI  
 ACCN: NM\_147025  
 Insert Size: 936 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_147025.1</a> , <a href="#">NP_667236.1</a>
<b>RefSeq Size:</b>	936 bp
<b>RefSeq ORF:</b>	936 bp
<b>Locus ID:</b>	259027
<b>Cytogenetics:</b>	11 B4
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]