

## Product datasheet for MC214187

### Olfr166 (NM\_147068) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Olfr166 (NM\_147068) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Olfr166  
**Synonyms:** MOR270-1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214187 representing NM\_147068  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGAGAAATGGAATCAGAGCTCAAGTATTTACTCTGTAGGACTGCTCCACAAAACCAACAGGCC  
 TGCTACTTTTATGCTCATCATCTTTGTCTTCTCTGGCTTTGTGTGGCAACTCAGGAATGATCCACCT  
 CATTTCGTGGATCCAAGGCTCCACACCCCATGACTTTCTCCTCAGTCAGCTCTCTCATGGACCTG  
 ATGTACATTTCTACCACTGTTCCCAAGATGGCATTAACTTCCTTTCTGGCCAGAAAAGCATCTCTTTTC  
 TGGGCTGTGGAGTGCAATCCTTCTTCTCCTGACTATGGCATGTTCTGAGGGCTTCTCTTGGCTCCAT  
 GGCTTATGATCGTTTTGTGGCTATCTGCCATCCCTTCACTATCCCATTCGCATGAGCAAAATAATGTGT  
 CTGAAGATGATCATAGGATCCTGGATATTGGGCTCAATCAACTCTTAGCACATACCGTCTATGCCCTTC  
 ATATTCCTTACTGCCATTCTAGGTCCATTAACCATTTCTCTGTGATGTTCCAGCCATGTTGCCCTTGGC  
 CTGTATGGACACTTGGGTTTATGAGTACATGGTGTGTTGTGAGCACAAGCCTGTTTCTCCTACTGCCTTC  
 CTTGGTATCACAGCTTCTATGGTCGGGTCCTTTTGTGTCTCCACATGCGCTCAAAAGAGGAAAGA  
 AGAAGGCCTTACCACATGCTCAACTCACTTAACTGTGGTGACATTTTACTATGCACCTTTTGTCTATAC  
 CTATCTTCGACCTAGGAGTCTTCGCTCCCAACAGAAGATAAGATTCTGGCTGTTTTCTACACTATCCTT  
 ACCCCCATGCTCAACCCATCATTTATAGTCTGAGGAATAAGGAGTCTGGGGCCATGACAAGAGTCC  
 TTGGTACTTTTCTTCAACTAAACCGTAA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_147068  
**Insert Size:** 939 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_147068.1</a> , <a href="#">NP_667279.1</a>
<b>RefSeq Size:</b>	939 bp
<b>RefSeq ORF:</b>	939 bp
<b>Locus ID:</b>	259071
<b>Cytogenetics:</b>	16 A3
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]