

Product datasheet for MC214297

Zfp446 (NM_001168561) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp446 (NM_001168561) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Zfp446
Synonyms:	A630035I11Rik; A930031C22; AA120513; AV306632
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC214297 representing NM_001168561 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCATCTCCTCGGGTACCCACACCTATCCCTCAGAGACTCCAGGACCACCCTTGAGGAGCCTGAGG
 CTGCACGACTCCGTTTCCGGGTTTCTGCTATGAGGAGGTTGAGGGACCCGAGAGGCACTGGCCAGCT
 TCGAGAGTTGTGTACCACTGGCTACAGCCTGAGTCATCCTCAAAGAACACATAATAGAGTTGCTGATT
 TTGGAACAGTTCCTGGGTGTGTGCCCCCTGAGATCCAGGCCTGGGTCGAGGGCAGCGACCAGGAAGCC
 CCGAGGAGGCCGAGTTCCTGTTGAGGGCTGCAGCATGATCCTGGGCAGCTGCTGGGCTGGGATCCGA
 GGGACATCTTGAACATCAGGAAACGGTCTCCATATCCTTCAAACAGGAAGGATCCAGGAGTGGGGCCTG
 TTGGAATCATCAAAAAAGAGCTGTGCTGGGTGTGATGCCGGAGAAGTACGACACAGTGGTGTCCAGG
 CAAGCCTGCCGCTGCTACAGCCAGAGACACATGTTGACTCTGAGCTGAGACCAAAACAGGAGATGCCACA
 CGAAGGACCAGAAATCCCTCCGGAGCCATCCCCCAGAAGTGGGAGCCATTGCAGATCCAAGGCTAGTTCAA
 GCCACTCCAAGTGAAAGGCAGAGCCCTTGAAGGACCCTTCAGTCTTGTCTCCACACCACTGCTTGAAG
 CCCCAGCGAGGCCACACCTCGGAAGCTGTACGTTTGTGAGCAGTGTGGCCTCAGCTTTGACTGGAAGTC
 CGTTTTTCATCATCACCTCCGCACACACGAGAGGGCCAGGCCTGGAGAGACCATCACAGGTGGCCCGG
 GAGCCGGCCATGAAGCGCTCCCTGGCCTTCGGGGATACACATGCGTGGAGTGTGGGCGCAGCTTCAGT
 GGAAGTCACAGCTGGTCAATCCACCAGAGCCATGCAGGCCAGCGGCCCACTTCTGCCCGACTGTGG
 CTGTAGCTTTGACTGGAAGTTTCAGCTGGTCATTACAGGAAGATCCACCAGCCAGAGGGTCCAT**GA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001168561
Insert Size:	1047 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001168561.1, NP_001162033.1</u>
RefSeq Size:	4498 bp
RefSeq ORF:	1047 bp
Locus ID:	269870
Cytogenetics:	7 A1