

Product datasheet for **MC214310**

Fam26d (NM_001081165) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam26d (NM_001081165) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fam26d
Synonyms: 4732454E20Rik; AI586164; Calhm4; Gm643
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC214310 representing NM_001081165
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCCCAGATCTCAACTGTATTTTCATCTTCTCTGCTCAGAAGTGAACCGTGCATCAATTCCTTTGATAG
CGATTTTGACTGTTTGGGACAGCAACTGTTCTCCTTATAACATTCAGTTGCCGTCAAGTCGGAAA
AAATTTCTATTATGGTTCAGCTTTTCTTGTGGTTCCTGCCTTGATCCTTCTGATTGCTGGCTATGCTCTA
CGAGGCCAAATGTGGACAGTTGCCAGCGAGTACTGCTGCTGCAGCTGTACCCCTCCATACCGGAGAAGCA
GCCCTCTGGAGAGGAGATTGGCTTGCCCTCATGTTCTTCGACATCACTGGGAGGGCTCTGGTTGCTCCACT
GACATGGTTGACGGTGACCCTGCTGACAGGTACCTACTATGAATGTGCAGCAAGTGAGTTTGCCTCTGTG
GACCAATACCCGATGTTTGCTAATGTCACCTCCTAGCAAACGAGAAGAGATGCTAGCTGGATTTCCATGCT
ACACGTCAGCTCCCTCTGATGTGATTCCAATAAGAGATGAAGTGGCTCTTCTACACAGATACCAAGTCA
AATGCTGGGCTGGATTCTGGTCGTTTTGGCAACCATGCTCTCCTGCTCTCCAAATGCTCTGGCAAGATGT
TGCTCTCCCTCACCTCTCTGCAGCATCACTATTGGACCAATCACCTCCACAACGAGAGAGTGTCTTTG
AGAAAGCTGCTGAGGAGCATTCCCAACTCTTATTAGACATCGAATCAAGAAAGTATTTGGCTTCGTTCC
TGGGAGTGAAGACATCAAACATATTCGATCCCTTCGTGTCAGGACTGGAGAGAAATCTCTGTACCCAAT
CTTCTATGTGTGGGTGACACCTCACAGGGTCCCTATAGCTTCTTGGAGACAGGGTGGTTGAGGAAAATG
AGGAAGATAGACAAGAAGGTATTGAAATGAAACCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2034_a12.zip

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001081165
Insert Size:	948 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001081165.1 , NP_001074634.1
RefSeq Size:	3714 bp
RefSeq ORF:	948 bp
Locus ID:	270711
UniProt ID:	Q8CE93
Cytogenetics:	10 B1
Gene Summary:	Pore-forming subunit of a voltage-gated ion channel.[UniProtKB/Swiss-Prot Function]