

Product datasheet for **MC214328**

Timd4 (NM_178759) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Timd4 (NM_178759) Mouse Untagged Clone
Tag: Tag Free
Symbol: Timd4
Synonyms: B430010N18Rik; TIM-4; Tim4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC214328 representing NM_178759
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCCAAGGGGCTTCTCCTCCTCTGGCTGGTGACGGAGCTCTGGTGGCTTTATCTGACACCAGCTGCCT
 CAGAGGATACAATAATAGGGTTTTGGCCAGCCGGTACTTTGCCTTGCATTACCTCTCGTGGTCCCA
 GAGCCGCAACAGTATGTGCTGGGCAAAGTTTCATGTCCCAATCCAAGTCAATGCAGAGCTTCTCCGT
 ACAGATGGAACAAGAATCATCTCCAGGAAGTCAACAAAATATACACTTTTGGGGAAGTCCAGTTTGGTG
 AAGTGTCTTGACCATCTCAAACACCAATCGAGGTGACAGTGGGGTGTACTGCTGCCGTATAGAGGTGCC
 TGGCTGGTTCAATGATGTCAAGAAGAATGTGCGCTTGGAGCTGAGGAGAGCCACAACAACCAAAAAACCA
 ACAACAACCACCCGCAACCAACCCCTTATGTGACCACCACCCAGAGCTGCTTCCAACAACAG
 TCATGACCACATCTGTTCTCCAACCAACCAACCCAGACTAGCCACCACTGCCTTCAGTACAGC
 AGTGACCACGTGCCCTCAACAACCTGGCTCCTTCTCAAGAACCACAAAAGGGTCCGCCTTCACT
 ACAGAATCAGAACTCTGCCTGCATCCAATCACTCTCAAAGAAGCATGATGACCATATCTACAGACATAG
 CCGTACTCAGGCCACAGGCTCTAACCTGGGATTCTCCCATCCACTTACAGCTGACGACACAGAAAAC
 AACATTAACAACAAGTGAAGTCTTTGCAGAAGACAATAATCACATCAGATCAACAGCAGACAGACCATC
 TTGATCATTGCCCTGCTGTGGGATTTGTGCTAATGGTGTATTGTTTCTGGCGTTTCTCCTTCGAGGGA
 AAGTCACAGGAGCCAAGTGTTCAGAGACACAAGAGGCCAGACAACACTGAAGATAGTGACAGCGTCCT
 CAATGACATGTACACGGGAGGGATGATGAAGACGGGATCTTCACTCTC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja3146_e05.zip



[View online »](#)

Restriction Sites: Sgfl-Mlul

ACCN: NM_178759

Insert Size: 1032 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178759.4](#), [NP_848874.3](#)

RefSeq Size: 2175 bp

RefSeq ORF: 1032 bp

Locus ID: 276891

UniProt ID: [Q6U7R4](#)

Cytogenetics: 11 B1.1

Gene Summary: Phosphatidylserine receptor that enhances the engulfment of apoptotic cells. Involved in regulating T-cell proliferation and lymphotoxin signaling. Ligand for HAVCR1/TIMD1. [UniProtKB/Swiss-Prot Function]