

Product datasheet for MC214337

Lkaaar1 (NM_199023) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Lkaaar1 (NM_199023) Mouse Untagged Clone
Tag: Tag Free
Symbol: Lkaaar1
Synonyms: 4930526D03Rik; Gm692
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214337 representing NM_199023
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGCCGACCTGGGAGTAAAGGGTGCACGGGAGCGGGACAAAAACAGTGCATCTGGAGCTGGAGCTGGAG
 CTGGAGCTGGAGCTGGAGCTGGAGAGAAGCACAGGAAGGGCCGAGGACTACAGACCCTCGAAGACAGG
 CTGGGCATTGACTAAGCAGAGGCTGGTAGCCTGTCCCCTACACTACGCCAACGCCACCTGCTCTTTGGC
 GACTTTCTTGACGACATAGGCAAGGTGGCTCCATGTTTCCACGCGAGTCGGTGGAGTTACCTACGACA
 TGCCCTGACCCGCGCACCTGGTACAGGCATTAATCTACCATCTGAGCATCAAAATCGGTTTCTGGGCC
 TATTAAGCCGCGGAAGCCCGGGCCGAGTCCATACCCTGCGGTTGCGCTACACCCGAATGCGGGCAGAG
 GAAATCTCTTTGCTGATCCAGAAGCAGAGCTCCGCGCGGCTGCCATCCGACTGGAATGTTCTACCAC
 CCCAACTGAAGCCCAAAAGATCCCTGATCCTCTGGACCGACACGAGCGGAGGCGTGTGAGACAATCCT
 GGAAGAGGAGGTAGATGGCAACATCTTCCCCGATGA

ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_199023
Insert Size: 597 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_199023.3</u> , <u>NP_950188.2</u>
RefSeq Size:	716 bp
RefSeq ORF:	597 bp
Locus ID:	277496
UniProt ID:	<u>Q8BIG2</u>
Cytogenetics:	2 H4