

Product datasheet for MC214340

Srarp (NM_001033374) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Srarp (NM_001033374) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Srarp
Synonyms:	Gm694
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC214340 representing NM_001033374 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTTTCTCAAAGGACCCAGGAGAACCAGCCTTAGGGATTCTAGTGTGGAGATGAGCTCTGGTACGC
 AGCCGTCTCGCACCCAAGGCCATGCCACAGCTCACGTCACTTTCCTTATTGACTGTGCAACTGGGAA
 ACAAGTCTCCCTGGCAGCATCCACGGCGCCACCCACGCATCAAGAGCCAATCAGGGATGTGTCGCCCCA
 CCCATGAAGACCTTCGTTATGTTCCGTGGAACCAACACGATGCTTGGGACTCAGAACATTTCTCTCAGTC
 GCGGAGCCCTTGACGGAGCCAAGGATACCTTGCCACCCTACAGAGGACTCGGGGCTCCACATTGCTCCC
 TGCCAGCCTACCTGGTCCCAGAATGATCCTAAAGCCCAGGGGAGCTCCCTGAAGCCTGGAGCTACTGAA
 AAGCACTCAACCCGGGAGAAGGTCAAACACTCACTGAAGGCTCTGACCTGTCTCTGCGGGCAGGTGGAGT
 AG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001033374
Insert Size:	492 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001033374.3 , NP_001028546.1
RefSeq Size:	922 bp
RefSeq ORF:	492 bp
Locus ID:	277744
UniProt ID:	Q3ULG3
Cytogenetics:	4 D3
Gene Summary:	May regulate the transcriptional function of androgen and estrogen receptors. [UniProtKB/Swiss-Prot Function]