

## Product datasheet for MC214377

### H2ac12 (NM\_175659) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** H2ac12 (NM\_175659) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** H2ac12  
**Synonyms:** Hist1h; Hist1h2ah  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Fully Sequenced ORF:** >MC214377 representing NM\_175659  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTCTGGACGCGGTAAGCAGGGCGGCAAGGCTCGCGCCAAGGCCAAGACCCGCTCCTCCGGGCCGGCC  
TGCAGTTCCTCCGTGGGCCGCGTGCACCGGCTGCTCCGCAAGGGCAACTACTCGGAGCGCTGGGCGCCGG  
CGCCCCGTGTACCTGGCGGCCGTGCTGGAGTACCTGACGGCCGAGATCCTGGAGCTGGCGGCAACGCG  
GCCCGGACAACAAGAAGACGCGCATCATCCCGGCCACCTGCAGCTGGCCATCCGCAACGACGAGGAGC  
TCAACAAGCTGCTGGGCCGCGTGACCATCGCGCAGGGCGGCGTCTGCCAACATCCAGGCCGTGCTGCT  
GCCAAGAAGACCGAGAGCCACCACAAGGCCAAG**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_175659  
**Insert Size:** 387 bp  
**OTI Disclaimer:**

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_175659.2</a> , <a href="#">NP_783590.1</a>
<b>RefSeq Size:</b>	522 bp
<b>RefSeq ORF:</b>	387 bp
<b>Locus ID:</b>	319168
<b>UniProt ID:</b>	<a href="#">Q8CGP6</a>
<b>Cytogenetics:</b>	13 A3.1
<b>Gene Summary:</b>	Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H2A family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. [provided by RefSeq, Aug 2015]