

Product datasheet for **MC214383**

H2bc1 (NM_175663) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: H2bc1 (NM_175663) Mouse Untagged Clone
Tag: Tag Free
Symbol: H2bc1
Synonyms: Hist1h; Hist1h2ba; Th2b
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC214383 representing NM_175663
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCGGAGGTGGCGGTAAGGGTACTACTTTCCAAGAAAGGCTTCAAGAAAGCGGTACCAAGACCC
AGAAAAAGGAGGGCCGAAACGTAAGAGATGCCGCAAGGAGAGCTACTCCATTTACATCTATAAGGTGCT
GAAACAAGTGCACCCCGACCCGGCATCTCTCCAAGGCCATGAGCATCATGAACCTTTGTGACAGAC
ATCTTCGAGCGCATCGCGAGCGAGGCGTCCCGCTGGCGCATTACAACAAGCGCTCGACCATCACGTCCC
GGGAGATCCAGACGGCCGTGCGCCTGCTGCTGCCCGGGAGCTGGCCAAGCACCGGTGTCGGAGGGCAC
CAAGGCCGTCACCAAGTACACCAGCTCCAAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_175663
Insert Size: 384 bp



[View online »](#)

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_175663.2](#), [NP_783594.1](#)

RefSeq Size: 432 bp

RefSeq ORF: 384 bp

Locus ID: 319177

UniProt ID: [P70696](#)

Cytogenetics: 13 A3.1

Gene Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H2A family. Transcripts from this gene contain a palindromic termination element. [provided by RefSeq, Aug 2015]